




Idys™-PLIF



Designed for the spine



CLARIANCE
innovation for spine surgery

Contents

Introduction	3
Features & benefits	4
Surgical technique	6
Product catalogue	14

Introduction

The Idys™-PLIF cage has been designed by a team with extensive experience in the development of spinal instrumentation. Made of PEEK, the Idys™-PLIF cage is specially designed for lumbar and lumbo-sacral interbody fusion via a bilateral posterior approach. The design team focused on creating a PLIF cage suited to the vertebral endplates anatomy that would maintain disc height, restore lumbar lordosis and stabilize the treated segment. Thus, Idys™-PLIF creates the most favorable conditions to optimize bone fusion. Its complete and ergonomic instrumentation meets the spine surgeon's need for reliability, safety and ease of use during surgical practice.

- » *Indications:*
- Lumbar degenerative discopathies and instabilities
 - Grade I or II spondylolisthesis, with or without stenosis
 - Pseudarthrosis or fusion failure

Caution : The use of the Idys™-PLIF cage must be systematically associated with a posterior fixation device (Erisma™-LP). For greater convenience, the pedicle screws can be inserted before or after the Idys™-PLIF cages are implanted.

Features & Benefits

Anatomical shape

▶ The convex shape of the Idys™-PLIF cage is perfectly adapted to the anatomy of the vertebral endplates. This shape allows restoration of the chosen disc height as well as lumbar lordosis of the treated segment. Close contact with the vertebral endplates provides effective integration and excellent stability of the implant within the intervertebral space.



Self-distracting tip

▶ The self-distracting tip of the Idys™-PLIF cage allows progressive distraction of the vertebral endplates in order to facilitate insertion and offer excellent penetration of the cage in the intervertebral space.



Grafting space maximized

▶ The size of the fusion chamber is optimized to offer maximal contact between vertebral endplates and grafts. Bone fusion is thereby promoted and fusion spaces are larger. Lateral side holes enable good graft vascularization.




X-ray markers (Tantalum)

▶ X-ray markers allow verification of the cage position during surgery and postoperative clinical follow-up. Fusion can be assessed and controlled during clinical follow-up thanks to the radiolucency of the PEEK material.



Ergonomic instrumentation

▶ The simple, intuitive and ergonomic instrumentation of the Idys™-PLIF cage guarantees accurate, quick and safe maneuvers to surgeons.



Versatile instrumentation

▶ Instrumentation is adaptable and can be used with different types of surgical techniques, including PLIF or PTLIF.



Surgical Technique


- | | |
|----|----------------------------------|
| 1 | Disc exposure |
| 2 | Distraction |
| 3 | Discectomy |
| 4 | Endplate preparation |
| 5 | Implant sizing |
| 6 | Implant preparation |
| 7 | Implant insertion |
| 8 | Insertion of the second cage |
| 9 | Final positioning |
| 10 | Supplementary posterior fixation |

In line with all necessary safety protocols, the patient is positioned on the operating table in the prone position. The final positioning of the patient and surgical approach are based on known techniques, routinely used by surgeons.

Disc exposure


1

Access to the disc space is obtained by performing a bilateral lamino-facetectomy using the **osteotome**. Resection must be sufficient to allow access of the instruments and insertion of the cage into the disc space.



It is advisable to preserve the facet joints in order to guarantee the stability of the treated segment and minimize disruption to adjacent segments.

The **nerve root retractor** is used to protect the surrounding nerve structures throughout the surgical procedure.



04781006
NERVE ROOT RETRACTOR

Distraction


2

Progressive distraction tension is applied in accordance with the surgeon's habits and preferences. This maneuver temporarily opens the posterior disc space and promotes increased exposure for discectomy, decompression and delivery of the implant.

Discectomy

3

A window is created in the intervertebral disc using a scalpel.




04705010 OSTEOTOME




04781006
NERVE ROOT RETRACTOR

Shavers are used in rotation to start the discectomy and gradually restore the disc height. Shavers come in sequential heights from 6mm to 14mm, increasing in 1 mm increments.




The disc material is removed using the **disc rongeur**.




When possible, the anterior and lateral walls of the annulus are preserved in order to provide additional stability for the two Idys™-PLIF cages.

Endplate preparation

The endplates are prepared using the **curette** and the **cup curette** to remove the remaining layers of the entire cartilaginous endplates and expose bleeding bone.



Thorough cleaning of the endplates is important for the vascular supply of the bone graft. The structural integrity of the endplates must be preserved to allow for structural support of the cage.



Implant sizing

The size and height of the final implant are determined by selecting the trial that restores the desired disc height and lumbar lordosis of the treated segment.

The **trial** is connected and secured to the **PLIF cage holder** by turning the locking knob clockwise.





04704015 CURETTE

04704005 CUP CURETTE

04714XXX PLIF TRIAL

04715002 PLIF CAGE HOLDER

The trial is then inserted in the intervertebral space using a surgical mallet. The trial is not dissociated from the PLIF cage inserter.



04762000 SLAP HAMMER

Distraction is momentarily released and an X-ray check is performed to verify the trial position and validate the size and height of the final implant.

Distraction is re-established and the trial is removed using the **slap hammer**.

The PLIF cage holder is disconnected from the trial by completely turning the locking knob.

Implant preparation

6

It is important to completely fill the implant cavity in order to ensure optimal contact between the graft material and the vertebral endplates.

The selected cage is connected to the **PLIF cage holder** and placed into the **PLIF filling block**.




The **PLIF graft impactor** is used to firmly pack the graft into the implant cavity.

Implant insertion

7

The prepared cage is inserted into the intervertebral space with light impactions. When inserting the cage, the **PLIF cage holder** must be turned medially to ensure that the tip of the cage is facing the entry point. During cage impaction, the PLIF cage holder must be held straight so that the implant is inserted into the disc space laterally.



04715002 PLIF CAGE HOLDER




04717003
PLIF FILLING BLOCK



0476100X
PLIF GRAFT IMPACTOR



An X-ray check is performed and the cage is disconnected.




Insertion of the second cage

8

The opposite side is prepared as described in steps 2 to 5 and the second cage is inserted similarly.



 Prior insertion of the second cage, graft can be inserted between the two implants.

Final positioning

9

If needed, the PLIF cage pusher may be used to adjust the position and the depth of the cage within the intervertebral space.





Distraction is released and a final X-ray check is performed.



Supplementary posterior fixation

10

Posterior fixation (Erisma™-LP) is required to place the treated segment under compression and to enhance the stability of the Idys™-PLIF cages.



04716001
PLIF CAGE PUSHER

Product Catalogue



OSTEOTOME 04705010



NERVE ROOT RETRACTOR 04781006



SHAVER
Ø 6mm 04707006
Ø 8mm 04707008
Ø 9mm 04707009
Ø 10mm 04707010
Ø 11mm 04707011
Ø 12mm 04707012
Ø 13mm 04707013
Ø 14mm 04707014



SLAP HAMMER 04762000



PLIF FILLING BLOCK 04717003



PLIF GRAFT IMPACTOR
L. 20mm 04761000
L. 25mm 04761001



STRAIGHT DISC RONGEUR 06708001



CURVED DISC RONGEUR 06708002



DISC RONGEUR 45° 06708003



PLIF CAGE PUSHER 04716001



CAGE PUSHER 04716002



CYLINDRICAL HANDLE
LARGE 99782003



CURETTE 04704015



CUP CURETTE 04704005




PLIF CAGE HOLDER 04715002




T-HANDLE 99781001


Product Catalogue




IDYS COMMON INSTRUMENTS TRAY 04990010



IDYS PLIF / PTLIF INSTRUMENTS TRAY 04990011
LID FOR IDYS PLIF / PTLIF INSTRUMENTS TRAY 04990013
PTLIF TRIALS TRAY 04990016
PLIF TRIALS TRAY 04990015





IDYS PLIF / PTLIF IMPLANTS TRAY 04990020
LID FOR IDYS PLIF / PTLIF IMPLANTS TRAY 04990027




TRIALS	L 20 MM
H 08 mm	04714208
H 09 mm	04714209
H 10 mm	04714210
H 11 mm	04714211
H 12 mm	04714212
H 13 mm	04714213

TRIALS	L 25 MM
H 08 mm	04714308
H 09 mm	04714309
H 10 mm	04714310
H 11 mm	04714311
H 12 mm	04714312
H 13 mm	04714313
H 14 mm	04714314



PLIF CAGE	L 20 MM
H 08 mm	04522008
H 09 mm	04522009
H 10 mm	04522010
H 11 mm	04522011
H 12 mm	04522012
H 13 mm	04522013

PLIF CAGE	L 25 MM
H 08 mm	04522508
H 09 mm	04522509
H 10 mm	04522510
H 11 mm	04522511
H 12 mm	04522512
H 13 mm	04522513
H 14 mm	04522514



The surgical technique shown is for illustrative purposes only. The technique actually employed in each case will always depend upon the medical judgment of the surgeon exercised before and during surgery as to the best mode of treatment for each patient. It is recommended to see the package insert for the complete list of indications, warnings, precautions, and other medical information.



CLARIANCE
18, Rue Robespierre
62217 Beaurains - FRANCE
Tel +33 (0)3 21 16 12 15
Fax +33 (0)3 21 15 50 73



18, Rue Robespierre
62217 Beaurains - FRANCE
[Tel +33 \(0\)3 21 16 12 15](tel:+33321161215)
[Fax +33 \(0\)3 21 15 50 73](tel:+33321155073)
contact@clariance-spine.com
www.clariance-spine.com

Distributed by