



STERILE REINFORCED SURGICAL GOWS

MDR 2017/745 Ref: 6200



Medical grade SMS

These extra protective gowns have an additional layer of waterproof material from the wrist to the elbow and on the front. Ideal for longer procedures and with large amounts of liquids. They are ideal for heavy procedures where surgeon contamination is crossed out and bidirectional and patient could occur.

CE Marked

STERILE | EO

Oley Sterile Reinforced Surgical Gowns are sterilized in accredited sterilization institute. Based on each shipment, we send sterilization report and our declaration.

SPECIFICATIONS

Disposable, sterile surgical gown for wet heavy surgical procedures

Made of non-woven fabric of the SMS type with a minimum weight of 75 gsm , reinforcement is made of made of PP / PE gr./m2

Comply with the standard EN 13795-1

Under the direction of the European union standard 93/42/CEE Class I

Expiration of the goods is in 5 years .

Manufactured in accordance with the normative requirements for medical devices

Wrapping cut, free neckline, covering the user's entire back

Fixation strap to secure the gown at the waist

Resistant to fluid leakage and resistant to tearing in both dry and wet conditions


Breathable and provides thermal physiological comfort

Long sleeves finished with knitted cuffs

Latex-free

Packed individually in a sterilization package

Gown fabric allows the skin to breathe, but it will remove alcohol, blood and similar body fluids from the environment and do not allow it to pass to the bottom



PACKING LIST

1 box has 50 pcs gowns

Box Dimension: 60*40*40 cm



Limited use



Do not dry



Do not wash



Do not dry-clean



Do not bleach



Do not iron



Flammable

Storage and disposal of accessories

- Store in dry and clean conditions, in the original packaging, away from direct sunlight, sources of high temperatures and solvent vapors
- Store at temperatures between -20 ° C and + 25 ° C (-4 ° F to + 68 ° F) and with relative humidity below 80%
- The expected shelf life is three years from the date of manufacture if stored as directed.
- Replace the product if it is damaged, heavily contaminated or in accordance with local regulations.
- Handle and dispose of contaminated product with care and in accordance with applicable regulations

