



varimed®

## PRODUCT CATALOGUE

PHLEBOLOGY · LYMPHOLOGY · ORTHOPAEDICS

VARİMED MEDİKAL SAN. VE TİC. LTD. ŞTİ.





## EDITORIAL

Varimed Medikal San. ve Tic. Ltd. Şti. was founded with the aim of medical supplies trade in 1999. Our company which always makes up a principle of customer satisfaction has decided to enter manufacturing of Compression Stockings, Anti-Embolism Stockings, Elastic Braces and Lymphoedema Armsleeves in 2008. For this purpose, considering recent technological developments, the installation of a modern plant was started to meet the requirements of Medical Device Directives and ISO 13485 Quality Management Systems in 3600 m<sup>2</sup> building in İstanbul-Bayrampaşa and our facility has got in to full swing in 2010. Our company with established manufacturing facility has come to position that can compete against the world's leading companies with it is quality and capacity. All stages from the raw material input to the product output are controlled with the excellent quality control system provide reaching users of the products as a flawless and try to offer a better product to the users that keep patient safety at the forefront with continuous R&D work. Today, Varimed has a wide range of products in phlebology, lymphology and orthopaedics.

Medical products of the Varimed meet all of the requirement of the European Directive 93/42/EEC and Turkish Medical Devices Directive. In accordance with regulations, products of Varimed that subject to the conformity assessment procedures are granted the CE mark.

Varimed compression stocking and garments, whether circular knit or flat knit, support stockings, anti-embolism stockings and supports regularly examined by independent test centres for any toxic substances according to the Oeko-Tex Standard 100. Granting of the Oeko-Tex Standard 100 Certificate guarantees tested quality for protection of your health.

Varimed maintains the Medical Devices Quality Management System according to DIN EN 13485. Since 2010, we demonstrated our ISO 13485 status to our customers with a certificate. This certification aims to ensure product traceability, effective production and risk management and compliance with regulations and customer requirements.

We aim to establish long-term customer relations based on trust and satisfaction with all our efforts.

**Dr. M. Volkan Öztürkatalay**  
CEO

# INDEX

## PHLEBOLOGY

<b>VARIMED MEDICAL COMPRESSION STOCKINGS .....</b>	<b>4-50</b>
Medical Compression Stocking Treatment in Venous Diseases .....	3-4
Varimed Premium .....	5-10
Varimed Classic .....	11-16
Varimed Cotton .....	17-22
Varimed Silver .....	23-28
Varimed Comfort.....	29-34
Varimed for Men .....	35-40
Varimed Strong .....	41-46
Varimed Leg Ulcer Kit .....	47-50

## **VARIMED SUPPORT STOCKINGS..... 51-64**

Support Stockings to Prevent Venous Diseases .....	51-52
Varimed Glamour 40 den.....	53-54
Varimed Glamour 70 den.....	55-56
Varimed Glamour 140 den.....	57-60
Varimed Style 280 den.....	61-64

## **VARIMED TRAVEL STOCKINGS ..... 65-68**

## **VARIMED ANTI-EMBOLISM STOCKINGS..... 69-71**

Varimed Classic Anti-embolism Stockings.....	72-74
Varimed Comfort Anti-embolism Stockings .....	75-78

## **VARIMED DIABETIC SOCKS ..... 79-82**



## **LYMPHOLOGY**

Compression Treatment in Lymphoedema .....	83-86
--	-------

## **LEG GARMENTS – ROUND KNIT ..... 87-92**

Varimed Strong Medical Compression Stocking .....	87-92
---	-------

## **LEG GARMENTS – FLAT KNIT ..... 93-102**

Varimed Forte Medical Compression Stocking .....	93-98
--	-------

Varimed Forte Toe Cap .....	99-102
-----------------------------	--------

## **ARM GARMENTS – ROUND KNIT ..... 103-118**

Varimed Strong Lymphoedema Armsleeves .....	103-106
---	---------

Varimed Strong Cotton Lymphoedema Armsleeves .....	107-112
--	---------

Varimed Strong Cotton Lymphoedema Gloves .....	113-118
--	---------

## **ARM GARMENTS – FLAT KNIT ..... 119-132**

Varimed Forte Lymphoedema Armsleeves.....	119-124
---	---------

Varimed Forte Lymphoedema Gloves.....	125-128
---------------------------------------	---------

Varimed Forte Seamless Lymphoedema Gloves .....	129-132
---	---------

## **ORTHOPAEDICS**

### **VARIMED ORTHOPAEDICS SUPPORTS ..... 133-161**

#### **BACK**

Varimed Lumbosacral Support .....	135-138
-----------------------------------	---------

Varimed Lumbosacral Support with Belt.....	135-138
--	---------

#### **ABDOMEN**

Varimed Abdominal Corset .....	139-142
--------------------------------	---------

## **KNEE..... 143-148**

Varimed Ligament and Patella Supported Knee Brace.....	144
Varimed Ligament and Patella Supported Knee Brace with Open Patella .....	144
Varimed Patella Supported Knee Brace.....	145
Varimed Knee Brace .....	145
Varimed Patella Strap.....	146

## **ANKLE..... 149-152**

Varimed Malleol Supported Ankle Brace.....	150
Varimed Ligament Supported Ankle Brace .....	150
Varimed Achilles Tendon Supported Ankle Brace.....	151
Varimed Ankle Brace .....	151

## **ARM ..... 153-154**

Varimed Shoulder Brace.....	154
-----------------------------	-----

## **ELBOW ..... 155-158**

Varimed Elbow Brace with silicone pads.....	156
Varimed Elbow Brace.....	156
Varimed Epycondylitis Bandage .....	157

## **HAND ..... 159-162**

Varimed Thumb Splint.....	160
Varimed Wrist Brace with Splint .....	160
Varimed Wrist Splint with Thumb Grip .....	161

## **VARIMED LIGHT SUPPORTS**

Sport Injuries and Elastic Braces.....	163-164
--	---------

## **VARIMED ELASTIC BRACES..... 165-168**

Varimed Elastic Ankle Brace .....	165
Varimed Elastic Calf Brace .....	165
Varimed Elastic Knee Brace .....	166
Varimed Elastic Elbow Brace.....	166



## **VARIMED DYNAMIC SUPPORTS ..... 169-176**

Varimed Dynamic Sport Stocking .....	170
Varimed Dynamic Back Support .....	171
Varimed Dynamic Knee Support .....	171
Varimed Dynamic Ankle Support .....	172
Varimed Dynamic Elbow Support .....	172
Varimed Dynamic Hand/Wrist Support .....	173
Varimed Dynamic Wrist Support.....	173

## **VARIMED THERMO SUPPORTS ..... 177-178**

Varimed Thermo Elastic Waist Corset .....	178
Varimed Thermo Knee Brace.....	178

## **ACCESSORIES ..... 179-183**

Varimed Donning Aid for Medical Compression Stocking .....	182
Varimed Dressing Aid for Medical Compression Stockings - Fabric.....	182
Varimed Dressing Aid for Compression Armsleeve - Fabric .....	183

## **SCIENTIFIC EVIDENCE OF COMPRESSION TREATMENT ..... 184-185**

## **CERTIFICATES ..... 186-189**

Oeko-Tex 100 Certificate .....	186
FDA Certificate .....	187
ISO 13485 Medical Devices Quality Management System Certificate .....	188
CE Certificate.....	189

## **MANUFACTURING LINE..... 190-192**







**PHLEBOLOGY**

## MEDICAL COMPRESSION STOCKING TREATMENT IN VEIN DISEASES



The veins return blood to the heart from all the organs of the body. The main problems of the veins include inflammation, clotting and defects that lead to distension and varicose veins.

The legs contain two major groups of veins: the superficial veins, located in the fatty layer under the skin and the deep veins, located in the muscles. Short veins connect the superficial and deep veins. Blood pressure in all veins is normally low and in the leg veins, low pressure can pose a problem. Blood has to flow from the leg veins upward to reach the heart when a person is standing. The deep veins play a major role in propelling blood upward. Located within the powerful calf muscles, these veins are forcefully compressed with every step.

Varicose veins are enlarged superficial veins in the legs. The precise cause of varicose veins is not known but it is probably a weakness in the walls of the superficial veins. This weakness may be inherited. Overtime, the weakness causes the veins to lose their elasticity.

When valves become incompetent the cycle of unidirectional blood flow is interrupted and backflow of blood occurs. This is most significant when the backflow blood occurs between the deep and superficial veins, as the increased pressure in the superficial veins will cause further valve incompetence. This is because the valve cusps no longer meet as a result of the stretching of veins. The overall effect of this increased superficial hydrostatic pressure is the formation of tortuous varicose veins.

If this increased hydrostatic pressure is not treated it leads to chronic venous insufficiency which in turn result in oedema, lipodermatosclerotic changes, varicose eczema and ulceration.

External compression applies controlled pressure to the skin. This then supports the superficial venous system, helping to reduce the increased hydrostatic pressure and redress the balance between the deep and superficial veins. This helps in reducing oedema and aiding venous return.

Unlike compression bandaging with its high working pressure, medical compression stockings exert a resting pressure. The daily build-up of pressure is controlled by the limited ability of the hosiery to stretch, so incompetent venous valves are approximated, venous return is accelerated, the fibrinolytic activity of the venous wall is increased and the risk of thrombosis reduced.

Medical Compression stockings exert an external pressure which is greater at ankle and reduces at the calf and thigh, thus increasing blood velocity within the deep venous system. How much pressure to apply is recognised that the amount of pressure required is dependent on the severity of the condition. Severe venous hypertension is associated with oedema, eczema, skin pigmentation, induration and ulceration. It can be seen that the management of these conditions will require a higher level of compression than conditions such as mild varicose veins and oedema.



## ■ Indications

### **Light (18- 21 mmHg) Compression Stockings:**

- For prevention of venous disease
- For slight varicosis without oedema and for early varicosis in pregnancy
- The diseases such as diabetes and rheumatism which lead to venous weakness and loss of elasticity
- For treatment after vein operations

### **Medium (23-32 mmHg) Compression Stockings:**

- For pronounced varicosis with a tendency to oedema
- For post-traumatic swelling
- For recovery period of simple ulcers
- For superficial thrombophlebitis
- For treatment after vein operations and sclerotherapy
- For pronounced varicosis in pregnancy

### **Strong (34-46 mmHg) Compression Stockings:**

- For severe conditional and post-thrombotic venous insufficiency
- For pronounced varicosis with oedema
- For secondary varicosis
- For dermatosclerosis
- For recovery period of serious and recurrent ulcers (after healing of ulcer)

### **Extra Strong (over 49 mmHg) Compression Stockings:**

- For lymphatic oedema
- For elephantiasis

## ■ Contraindications:

### **Absolute contraindications:**

- Progressive periphery arterial congestion disease
- Decompensated heart diseases
- Septic phlebitis
- Phlegmasia coerulea dolens

### **Relative contraindications:**

- Wetting dermatosis
- Incompatibility to compression stocking fabric
- Sensibility disturbances of the limb
- Progressive periphery neuropathy

VARIMED Premium



# ■ VARIMED Premium MEDICAL COMPRESSION STOCKINGS

Fashionable appearance

## Features:

Varimed Premium offer the ladies a compression stocking in which they feel beautiful all the time. This line brings together the requirements of medical compression stocking and fashion. You'll find a high therapeutic efficacy and elegant style, both in one product.

- CCL 1, 2
- Fashionable colours
- Slight transparent appearance
- Fineness
- Very soft and smooth touch
- Microfilament yarns give the stocking excellent breathability
- High moisture transportation for pleasant wearing comfort.
- Perfect fit because of soft stretch material. It has been manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Attractive topband
- Easy put on and take off
- Freedom of movements
- Soft, elastic toe and heel
- High wearing comfort
- Perfect sizing. Seven sizes
- For patients with slightly soft connective tissue
- Skin friendly - Tested to Oeko-Tex Standard 100
- Black, Dark Anthracite, Navy, Brown, Caramel, Burgund



LYCRA® elyaf INVISTA'nın ticari markasıdır.  
LYCRA® fibre is a registered trademark of INVISTA.

## ■ COLORS



**Black**



**Dark Anthracite**



**Navy**



**Brown**



**Caramel**



**Burgundy**





**Thigh Length Stocking With Fantasy Topband**



**Panty Hose Stocking**



**Panty Hose Stocking**



## ■ PRODUCT TYPES AND SIZES FOR VARIMED Premium

	STOCKING TYPE				FOOT		LENGTH		SIZE						
Compression Class	AD	AG	AT	ATU	OPEN	CLOSED	SHORT	NORMAL	1	2	3	4	5	6	7
I (18-21 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
II (23-32 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●

**AD:** KNEE LENGTH COMPRESSION STOCKINGS

**AG:** THIGH LENGTH COMPRESSION STOCKINGS

**AT:** PANTY HOSE COMPRESSION STOCKINGS

**ATU:** MATERNITY PANTY HOSE COMPRESSION STOCKINGS

● Available

○ Not Available

## ■ CODES AND DESCRIPTIONS OF VARIMED Premium

### Ccl I LIGHT COMPRESSION STOCKINGS (18 - 21 mmHg)


CODE	PRODUCT
3111	Light Compression Knee Length (AD) Closed Toe Stocking
3131	Light Compression Thigh Length (AG) Closed Toe Stocking
3141	Light Compression Panty Hose (AT) Closed Toe Stocking
3151	Light Compression Maternity Panty Hose (ATU) Closed Toe Stocking
3110	Light Compression Knee Length (AD) Open Toe Stocking
3130	Light Compression Thigh Length (AG) Open Toe Stocking
3140	Light Compression Panty Hose (AT) Open Toe Stocking
3150	Light Compression Maternity Panty House (ATU) Open Toe Stocking

### Ccl II MEDIUM COMPRESSION STOCKINGS (23 - 32 mmHg)

CODE	PRODUCT
3211	Medium Compression Knee Length (AD) Closed Toe Stocking
3231	Medium Compression Thigh Length (AG ) Closed Toe Stocking
3241	Medium Compression Panty Hose (AT) Closed Toe Stocking
3251	Medium Compression Maternity Panty Hose (ATU) Closed Toe Stocking
3210	Medium Compression Knee Length (AD) Open Toe Stocking
3230	Medium Compression Thigh Length (AG) Open Toe Stocking
3240	Medium Compression Panty Hose (AT) Open Toe Stocking
3250	Medium Compression Maternity Panty Hose (ATU) Open Toe Stocking

## ■ SIZE TABLE

	SIZE						
Leg Circumference (cm)	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cB1	23-27	24-29	26-32	29-35	31-37	33-39	35-41
cC	28-34	30-37	33-40	35-43	37-45	39-47	41-49
cD	27-33	29-36	32-39	34-42	36-44	38-46	40-48
cE	32-38	34-41	37-44	39-47	41-49	43-51	45-53
cF	37-47	40-52	43-56	46-60	50-65	53-68	56-71
cG	43-57	46-62	49-67	52-72	56-76	60-79	64-82



Length of Leg (cm)	
lG Normal	71-82
lG Short	62-70
lD Normal	38-43
lD Short	33-37

**Choosing the Size:** Two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg (lD is for knee-length stockings, lF is for mid-thigh length stockings, lG is for thigh-length and pantyhose stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.



## VARIMED Classic



# ■ VARİMED classic MEDICAL COMPRESSION STOCKINGS

## Elegant and Robust Compression Stocking

Varimed Classic Medical Compression Stockings are made from high quality microfilament yarn and medical type Lycra® which has a soft fabric that breathes on the skin and offer perfect fit and excellent wearing comfort. Due to the firm knitted fabric this quality is optimal for patients with soft connective tissue and for the treatment of serious indications.

### Features:

- CCL 1, 2, 3, 4
- Elegant appearance
- Microfilament yarns give the stocking excellent breathability
- High moisture transportation for pleasant wearing comfort.
- Perfect fit because of soft stretch material. It has been manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Easy put on and take off
- Freedom of movements
- Soft, elastic toe and heel
- High wearing comfort and high efficiency
- Perfect sizing. Seven sizes and two length
- Soft and smooth touch
- Opaque appearance
- For patients with soft connective tissue
- For patient with more severe phlebological indications
- Skin friendly - Tested to Oeko-Tex Standard 100




LYCRA® elyaf INVISTA'nın ticari markasıdır.  
LYCRA® fibre is a registered trademark of INVISTA.

## ■ COLORS



**Black**



**Beige**



***Knee Length Stocking***



***Thigh Length Stocking***



***Open Toe***



***Panty Hose Stocking***



## ■ PRODUCT TYPES AND SIZES FOR VARIMED Classic

	STOCKING TYPE					FOOT		LENGTH			SIZE								
Compression Class	AD	AG	AG/T	AT	ATU	OPEN	CLOSED	SHORT	NORMAL	LONG	0	1	2	3	4	5	6	7	8
I (18-21 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
II (23-32 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
III (34-46 mmHg)	●	●	●	●	○	●	○	●	●	●	●	●	●	●	●	●	●	●	●
IV (over 49mmHg)	●	●	○	○	○	●	○	●	●	●	○	●	●	●	●	●	●	●	○

**AD:** KNEE LENGTH COMPRESSION STOCKINGS

**AG:** THIGH LENGTH COMPRESSION STOCKINGS

**AG/T:** THIGH LENGTH WITH WAIST BELT COMPRESSION STOCKINGS

**AT:** PANTY HOSE COMPRESSION STOCKINGS

**ATU:** MATERNITY PANTY HOSE COMPRESSION STOCKINGS

● Available

○ Not Available

## ■ CODES AND DESCRIPTIONS OF VARIMED Classic

### Ccl I LIGHT COMPRESSION STOCKINGS (18-21 mmHg)

CODE	PRODUCT
<b>9111</b>	Light Compression Knee Length (AD) Closed Toe Stocking
<b>9131</b>	Light Compression Thigh Length (AG) Closed Toe Stocking
<b>91R1</b>	Light Compression Thigh Length with Waist Belt Right Leg (AG/T) Closed Toe Stocking
<b>91L1</b>	Light Compression Thigh Length with Waist Belt Left Leg (AG/T) Closed Toe Stocking
<b>9141</b>	Light Compression Panty Hose (AT) Closed Toe Stocking
<b>9151</b>	Light Compression Maternity Panty Hose (ATU) Closed Toe Stocking
<b>9110</b>	Light Compression Knee Length (AD) Open Toe Stocking
<b>9130</b>	Light Compression Thigh Length (AG) Open Toe Stocking
<b>91R0</b>	Light Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
<b>91L0</b>	Light Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
<b>9140</b>	Light Compression Panty Hose (AT) Open Toe Stocking
<b>9150</b>	Light Compression Maternity Panty Hose (ATU) Open Toe Stocking



## Ccl II MEDIUM COMPRESSION STOCKINGS (23-32 mmHg)

CODE	PRODUCT
9211	Medium Compression Knee Length (AD) Closed Toe Stocking
9231	Medium Compression Thigh Length (AG) Closed Toe Stocking
92R1	Medium Compression Thigh Length with Waist Belt Right Leg (AG/T) Closed Toe Stocking
92L1	Medium Compression Thigh Length with Waist Belt Left Leg (AG/T) Closed Toe Stocking
9241	Medium Compression Panty Hose (AT) Closed Toe Stocking
9251	Medium Compression Maternity Panty Hose Closed Toe (ATU) Stocking
9210	Medium Compression Knee Length (AD) Open Toe Stocking
9230	Medium Compression Thigh Length (AG) Open Toe Stocking
92R0	Medium Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
92L0	Medium Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
9240	Medium Compression Panty Hose (AT) Open Toe Stocking
9250	Medium Compression Maternity Panty Hose (ATU) Open Toe Stocking

## Ccl III STRONG COMPRESSION STOCKINGS (34 -46 mmHg)


CODE	PRODUCT
9310	Strong Compression Knee Length (AD) Open Toe Stocking
9330	Strong Compression Thigh Length (AG) Open Toe Stocking
93R0	Strong Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
93L0	Strong Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
9340	Strong Compression Panty Hose (AT) Open Toe Stocking

## Ccl IV EXTRA STRONG COMPRESSION STOCKINGS (over 49 mmHg )

CODE	PRODUCT
9410	Extra strong Compression Knee Length (AD) Open Toe Stocking
9430	Extra strong Compression Thigh Length (AG) Open Toe Stocking

## ■ SIZE TABLE

	SIZE								
Leg Circumference (cm)	0	1	2	3	4	5	6	7	8
cB	16-18	18-20	20-22	22-24	24-26	26-28	28-30	30-32	32-35
cB1	21-25	23-27	24-29	26-32	29-35	31-37	33-39	35-41	37-44
cC	26-31	28-34	30-37	33-40	35-43	37-45	39-47	41-49	43-52
cD	25-30	27-33	29-36	32-39	34-42	36-44	38-46	40-48	42-51
cE	29-35	32-38	34-41	37-44	39-47	41-49	43-51	45-53	47-55
cF	34-43	37-47	40-52	43-56	46-60	50-65	53-68	56-71	60-78
cG	40-53	43-57	46-62	49-67	52-72	56-76	60-79	64-82	69-92



Length of Leg (cm)	
lG Normal	71-82
lG Short	62-70
lD Normal	38-43
lD Short	33-37

**Choosing the Size:** Two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg (lD is for knee-length stockings, lF is for mid-thigh length stockings, lG is for thigh-length and pantyhose stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.

## VARIMED Cotton





# ■ VARIMED Cotton MEDICAL COMPRESSION STOCKINGS

## Perfect solution for sensitive skin

Varimed Cotton Compression Socks are comfortably worn by patients with sensitive skin or allergy sufferers due to skin compatibility. The inner surface of the sock that touches the skin is 100% cotton. Due to the firm knitted fabric this quality is optimal for patients with soft connective tissue and for the treatment of serious indications.

## Features:

- CCL 1, 2, 3
- Elegant appearance
- Extremely robust
- Microfilament yarns give the stocking excellent breathability
- High moisture transportation for pleasant wearing comfort
- Perfect fit because of soft stretch material. It has been manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Easy put on and take off
- Freedom of movements
- Soft, elastic toe and heel
- High wearing comfort and high efficiency
- Perfect sizing. Seven sizes and two length
- Soft and smooth touch
- Opaque appearance
- Particularly suitable for patients with soft connective tissue
- For patient with more severe phlebological indications
- For patient with mild lymphological indications
- Skin friendly - Tested to Oeko-Tex Standard 100



LYCRA® elyaf INVISTA'nın ticari markasıdır.  
LYCRA® fibre is a registered trademark of INVISTA.

## ■ COLORS



**Black**



**Beige**



***Knee Length Stocking***



***Thigh Length Stocking  
With Dot Silicone Topband***



***Thigh Length Stocking***

## ■ PRODUCT TYPES AND SIZES FOR VARIMED Cotton

Compression Class	STOCKING TYPE					FOOT		LENGTH			SIZE						
	AD	AG	AG/T	AT	ATU	OPEN	CLOSED	SHORT	NORMAL	LONG	1	2	3	4	5	6	7
I (18-21 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
II (23-32 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
III (34-46 mmHg)	●	●	●	●	○	●	○	●	●	●	●	●	●	●	●	●	●

**AD:** KNEE LENGTH COMPRESSION STOCKINGS

**AG:** THIGH LENGTH COMPRESSION STOCKINGS

**AG/T:** THIGH LENGTH WITH WAIST BELT COMPRESSION STOCKINGS

**AT:** PANTY HOSE COMPRESSION STOCKINGS

**ATU:** MATERNITY PANTY HOSE COMPRESSION STOCKINGS

● Available

○ Not Available

## ■ CODES AND DESCRIPTIONS OF VARIMED Cotton

### Ccl I LIGHT COMPRESSION STOCKINGS (18 – 21 mmHg)

CODE	PRODUCT
<b>8111</b>	Light Compression Knee Length (AD) Closed Toe Stocking
<b>8131</b>	Light Compression Thigh Length (AG) Closed Toe Stocking
<b>81R1</b>	Light Compression Thigh Length with Waist Belt Right Leg (AG/T) Closed Toe Stocking
<b>81L1</b>	Light Compression Thigh Length with Waist Belt Left Leg (AG/T) Closed Toe Stocking
<b>8141</b>	Light Compression Panty Hose (AT) Closed Toe Stocking
<b>8151</b>	Light Compression Maternity Panty Hose (ATU) Closed Toe Stocking
<b>8110</b>	Light Compression Knee Length (AD) Open Toe Stocking
<b>8130</b>	Light Compression Thigh Length (AG) Open Toe Stocking
<b>81R0</b>	Light Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
<b>81L0</b>	Light Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
<b>8140</b>	Light Compression Panty Hose (AT) Open Toe Stocking
<b>8150</b>	Light Compression Maternity Panty Hose (ATU) Open Toe Stocking

## Ccl II MEDIUM COMPRESSION STOCKINGS (23 – 32 mmHg)

CODE	PRODUCT
<b>8211</b>	Medium Compression Knee Length (AD) Closed Toe Stocking
<b>8231</b>	Medium Compression Thigh Length (AG) Closed Toe Stocking
<b>82R1</b>	Medium Compression Thigh Length with Waist Belt Right Leg (AG/T) Closed Toe Stocking
<b>82L1</b>	Medium Compression Thigh Length with Waist Belt Left Leg (AG/T) Closed Toe Stocking
<b>8241</b>	Medium Compression Panty Hose (AT) Closed Toe Stocking
<b>8251</b>	Medium Compression Maternity Panty Hose (ATU) Closed Toe Stocking
<b>8210</b>	Medium Compression Knee Length (AD) Open Toe Stocking
<b>8230</b>	Medium Compression Thigh Length (AG) Open Toe Stocking
<b>82R0</b>	Medium Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
<b>82L0</b>	Medium Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
<b>8240</b>	Medium Compression Panty Hose (AT) Open Toe Stocking
<b>8250</b>	Medium Compression Maternity Panty Hose (ATU) Open Toe Stocking

## Ccl III STRONG COMPRESSION STOCKINGS (34 -46 mmHg)


CODE	PRODUCT
<b>8310</b>	Strong Compression Knee Length (AD) Open Toe Stocking
<b>8330</b>	Strong Compression Thigh Length (AG) Open Toe Stocking
<b>83R0</b>	Strong Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
<b>83L0</b>	Strong Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
<b>8340</b>	Strong Compression Panty Hose (AT) Open Toe Stocking





## ■ SIZE TABLE

	SIZE						
Leg Circumference (cm)	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cB1	23-27	24-29	26-32	29-35	31-37	33-39	35-41
cC	28-34	30-37	33-40	35-43	37-45	39-47	41-49
cD	27-33	29-36	32-39	34-42	36-44	38-46	40-48
cE	32-38	34-41	37-44	39-47	41-49	43-51	45-53
cF	37-47	40-52	43-56	46-60	50-65	53-68	56-71
cG	43-57	46-62	49-67	52-72	56-76	60-79	64-82



Length of Leg (cm)	
lG Normal	71-82
lG Short	62-70
lD Normal	38-43
lD Short	33-37

**Choosing the Size:** Two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg (lD is for knee-length stockings, lF is for mid-thigh length stockings, lG is for thigh-length and pantyhose stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.

VARIMED Silver



# ■ VARIMED Silver MEDICAL COMPRESSION STOCKINGS

## Antibacterial Compression stocking with cotton and silver yarn

### Features of X-static silver yarn:

- Antibacterial. Fights bacteria which cause odours and eliminates fungi responsible for feet infection. Promotes hygiene of the feet and legs
- Antistatic. Disperses the electrostatic charge
- Therapeutic. Thanks to its conductive characteristics this new fibre gives benefits to health and favours blood circulation. It is also good for diabetic patients.
- Thermodynamic. The reflective and radiating qualities of the material encourage thermal conductivity ensuring that the stockings remain fresh in summer and warm in winter.
- Anti-odour. Neutralizes ammonia and denatured proteins, two principle causes of odour.

### Features of Compression Stocking:

- CCL 1, 2, 3
- Elegant appearance
- Microfilament yarns give the stocking excellent breathability
- High moisture transportation for pleasant wearing comfort
- Perfect fit because of soft stretch material. It has been manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Easy put on and take off
- Freedom of movements
- Soft, elastic toe and heel
- High wearing comfort and high efficiency
- Perfect sizing. Seven sizes
- Soft and smooth touch
- Opaque appearance
- For patients with soft connective tissue
- For patient with more severe phlebological indications
- Skin friendly - Tested to Oeko-Tex Standard 100



**X-static®**  
THE SILVER FIBER™



LYCRA® elyaf INVISTA'nın ticari markasıdır.  
LYCRA® fibre is a registered trademark of INVISTA.

## ■ COLORS



**Black**

**Beige**



***Thigh Length Stocking***



***Panty Hose Stocking***



***Thigh Length Stocking***



## ■ PRODUCT TYPES AND SIZES FOR VARIMED Silver

Compression Class	STOCKING TYPE				FOOT		LENGTH		SIZE						
	AD	AG	AT	ATU	OPEN	CLOSED	SHORT	NORMAL	1	2	3	4	5	6	7
I (18 - 21 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
II (23 - 32 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
III (34 - 46 mmHg)	●	●	●	○	●	○	●	●	●	●	●	●	●	●	●

**AD:** KNEE LENGTH COMPRESSION STOCKINGS

**AG:** THIGH LENGTH COMPRESSION STOCKINGS

**AT:** PANTY HOSE COMPRESSION STOCKINGS

**ATU:** MATERNITY PANTY HOSE COMPRESSION STOCKINGS

● Available

○ Not Available

## ■ CODES AND DESCRIPTIONS OF VARIMED Silver

### Ccl I LIGHT COMPRESSION STOCKINGS (18 - 21 mmHg)

CODE	PRODUCT
<b>10111</b>	Light Compression Knee Length (AD) Closed Toe Stocking
<b>10131</b>	Light Compression Thigh Length (AG) Closed Toe Stocking
<b>10141</b>	Light Compression Panty Hose (AT) Closed Toe Stocking
<b>10151</b>	Light Compression Maternity Panty Hose (ATU) Closed Toe Stocking
<b>10110</b>	Light Compression Knee Length (AD) Open Toe Stocking
<b>10130</b>	Light Compression Thigh Length (AG) Open Toe Stocking
<b>10140</b>	Light Compression Panty Hose (AT) Open Toe Stocking
<b>10150</b>	Light Compression Maternity Panty Hose (ATU) Open Toe Stocking

### Ccl II MEDIUM COMPRESSION STOCKINGS (23 – 32 mmHg)


CODE	PRODUCT
<b>10211</b>	Medium Compression Knee Length (AD) Closed Toe Stocking
<b>10231</b>	Medium Compression Thigh Length (AG) Closed Toe Stocking
<b>10241</b>	Medium Compression Panty Hose (AT) Closed Toe Stocking
<b>10251</b>	Medium Compression Maternity Panty Hose (ATU) Closed Toe Stocking
<b>10210</b>	Medium Compression Knee Length (AD) Open Toe Stocking
<b>10230</b>	Medium Compression Thigh Length (AG) Open Toe Stocking
<b>10240</b>	Medium Compression Panty Hose (AT) Open Toe Stocking
<b>10250</b>	Medium Compression Maternity Panty Hose (ATU) Open Toe Stocking

### Ccl III STRONG COMPRESSION STOCKINGS (34 -46 mmHg)

CODE	PRODUCT
<b>10310</b>	Strong Compression Knee Length (AD) Open Toe Stocking
<b>10330</b>	Strong Compression Thigh Length (AG) Open Toe Stocking
<b>10340</b>	Strong Compression Panty Hose (AT) Open Toe Stocking

## ■ SIZE TABLE

	SIZE						
Leg Circumference (cm)	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cB1	23-27	24-29	26-32	29-35	31-37	33-39	35-41
cC	28-34	30-37	33-40	35-43	37-45	39-47	41-49
cD	27-33	29-36	32-39	34-42	36-44	38-46	40-48
cE	32-38	34-41	37-44	39-47	41-49	43-51	45-53
cF	37-47	40-52	43-56	46-60	50-65	53-68	56-71
cG	43-57	46-62	49-67	52-72	56-76	60-79	64-82



Length of Leg (cm)	
lG Normal	71-82
lG Short	62-70
lD Normal	38-43
lD Short	33-37

**Choosing the Size:** Two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg (lD is for knee-length stockings, lF is for mid-thigh length stockings, lG is for thigh-length and pantyhose stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.



VARIMED Comfort





# ■ VARIMED Comfort MEDICAL COMPRESSION STOCKINGS

## Economic and reliable compression treatment for daily use


Varimed Comfort Compression Stockings are suitable for everyday use due to its robust and skin friendly material. Through carefully selected materials and knitting technology used Varimed Comfort Line is at the same time soft and resistant.

### Features:

- CCL 1, 2, 3
- High wearing comfort and efficacy
- Soft touch
- Opaque appearance
- Extremely durable yarns
- Easy put on and take off
- Freedom of movements
- Soft, elastic toe and heel
- Perfect fit
- Perfect sizing. Seven sizes
- For patients with slightly soft connective tissue
- For patient with more severe phlebological indications
- Skin friendly - Tested to Oeko-Tex Standard 100



## ■ COLORS



*Closed Toe*



*Open Toe*



***Panty Hose Stocking***



***Knee Length Stocking***



***Maternity Panty Hose Stocking***

## ■ PRODUCT TYPES AND SIZES FOR VARIMED Comfort

Compression Class	STOCKING TYPE					FOOT		LENGTH		SIZE						
	AD	AG	AG/T	AT	ATU	OPEN	CLOSED	SHORT	NORMAL	1	2	3	4	5	6	7
I (18 - 21 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
II (23 - 32 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
III (34 - 46 mmHg)	●	●	●	●	○	●	○	●	●	●	●	●	●	●	●	●

**AD:** KNEE LENGTH COMPRESSION STOCKINGS

**AG:** THIGH LENGTH COMPRESSION STOCKINGS

**AG/T:** THIGH LENGTH WITH WAIST BELT COMPRESSION STOCKINGS

**AT:** PANTY HOSE COMPRESSION STOCKINGS

**ATU:** MATERNITY PANTY HOSE COMPRESSION STOCKINGS

● Available

○ Not Available

## ■ CODES AND DESCRIPTIONS OF VARIMED Comfort

### Ccl I LIGHT COMPRESSION STOCKINGS (18 – 21 mmHg)

CODE	PRODUCT
<b>7111</b>	Light Compression Knee Length (AD) Closed Toe Stocking
<b>7131</b>	Light Compression Thigh Length (AG) Closed Toe Stocking
<b>71R1</b>	Light Compression Thigh Length with Waist Belt Right Leg (AG/T) Closed Toe Stocking
<b>71L1</b>	Light Compression Thigh Length with Waist Belt Left Leg (AG/T) Closed Toe Stocking
<b>7141</b>	Light Compression Panty Hose (AT) Closed Toe Stocking
<b>7151</b>	Light Compression Maternity Panty Hose (ATU) Closed Toe Stocking
<b>7110</b>	Light Compression Knee Length (AD) Open Toe Stocking
<b>7130</b>	Light Compression Thigh Length (AG) Open Toe Stocking
<b>71R0</b>	Light Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
<b>71L0</b>	Light Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
<b>7140</b>	Light Compression Panty Hose (AT) Open Toe Stocking
<b>7150</b>	Light Compression Maternity Panty Hose (ATU) Open Toe Stocking

## Ccl II MEDIUM COMPRESSION STOCKINGS (23 – 32 mmHg)

CODE	PRODUCT
7211	Medium Compression Knee Length (AD) Closed Toe Stocking
7231	Medium Compression Thigh Length (AG) Closed Toe Stocking
72R1	Medium Compression Thigh Length with Waist Belt Right Leg(AT) Closed Toe Stocking
72L1	Medium Compression Thigh Length with Waist Belt Left Leg(AT) Closed Toe Stocking
7241	Medium Compression Panty Hose (AT) Closed Toe Stocking
7251	Medium Compression Maternity Panty Hose (ATU) Closed Toe Stocking
7210	Medium Compression Knee Length (AD) Open Toe Stocking
7230	Medium Compression Thigh Length (AG) Open Toe Stocking
72R0	Medium Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
72L0	Medium Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
7240	Medium Compression Panty Hose (AT) Open Toe Stocking
7250	Medium Compression Maternity Panty Hose (ATU) Open Toe Stocking

## Ccl III STRONG COMPRESSION STOCKINGS (34 -46 mmHg)

CODE	PRODUCT
7310	Strong Compression Knee Length (AD) Open Toe Stocking
7330	Strong Compression Thigh Length (AG) Open Toe Stocking
73R0	Strong Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
73L0	Strong Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
7340	Strong Compression Panty Hose (AT) Open Toe Stocking



**Thigh Length Stocking**




**Knee Length Stocking**



## ■ SIZE TABLE

	SIZE						
Leg Circumference (cm)	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cB1	23-27	24-29	26-32	29-35	31-37	33-39	35-41
cC	28-34	30-37	33-40	35-43	37-45	39-47	41-49
cD	27-33	29-36	32-39	34-42	36-44	38-46	40-48
cE	32-38	34-41	37-44	39-47	41-49	43-51	45-53
cF	37-47	40-52	43-56	46-60	50-65	53-68	56-71
cG	43-57	46-62	49-67	52-72	56-76	60-79	64-82



Length of Leg (cm)	
lG Normal	71-82
lG Short	62-70
lD Normal	38-43
lD Short	33-37

**Choosing the Size:** Two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg (lD is for knee-length stockings, lF is for mid-thigh length stockings, lG is for thigh-length and pantyhose stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.

## VARIMED For Men



## ■ VARIMED For Men MEDICAL COMPRESSION STOCKINGS

Medical Compression Stocking for men with high wearing comfort. Varimed for Men Compression Stockings are exactly the same characteristics with Classic and Cotton Series and in addition to that it has ribbed design specially developed for men.

- CCL 1,2,3
- Ribbed design men stocking under Classic line and Cotton Line
- Below knee and Thigh Length
- Microfilament yarns give the stocking excellent breathability
- High moisture transportation for pleasant wearing comfort.
- Perfect fit thanks to soft stretch material, manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Easy put on and take off
- Freedom of movements
- Soft, elastic toe and heel
- High wearing comfort and high efficiency
- Perfect sizing. 7 sizes
- Opaque appearance
- Skin friendly - Tested to Oeko-Tex Standard 100



LYCRA® elyaf INVISTA'nın ticari markasıdır.  
LYCRA® fibre is a registered trademark of INVISTA.

## ■ COLORS



**Black**

## ■ PRODUCT TYPES AND SIZES FOR VARIMED For Men

	STOCKING TYPE		FOOT		LENGTH		SIZE						
Compression Class	AD	AG	OPEN	CLOSED	SHORT	NORMAL	1	2	3	4	5	6	7
I (18 - 21 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●
II (23 - 32 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●
III (34 - 46 mmHg)	●	●	●	○	●	●	●	●	●	●	●	●	●

**AD:** KNEE LENGTH COMPRESSION STOCKINGS

**AG:** THIGH LENGTH COMPRESSION STOCKINGS

● Available

○ Not Available

## ■ CODES AND DESCRIPTIONS OF VARIMED For Men

### VARIMED Classic For Men

#### Ccl I LIGHT COMPRESSION STOCKINGS (18 – 21 mmHg)

CODE	PRODUCT
<b>9181</b>	Light Compression Knee Length (AD) Closed Toe Stocking
<b>9191</b>	Light Compression Thigh Length (AG) Closed Toe Stocking
<b>9180</b>	Light Compression Knee Length (AD) Open Toe Stocking
<b>9190</b>	Light Compression Thigh Length (AG) Open Toe Stocking

### VARIMED Cotton For Men

#### Ccl I LIGHT COMPRESSION STOCKINGS (18 – 21 mmHg)

CODE	PRODUCT
<b>8181</b>	Light Compression Knee Length (AD) Closed Toe Stocking with cotton
<b>8191</b>	Light Compression Thigh Length (AG) Closed Toe Stocking with cotton
<b>8180</b>	Light Compression Knee Length (AD) Open Toe Stocking with cotton
<b>8190</b>	Light Compression Thigh Length (AG) Open Toe Stocking with cotton



## VARIMED Classic For Men

### Ccl II MEDIUM COMPRESSION STOCKINGS (23 – 32 mmHg)

CODE	PRODUCT
9281	Medium Compression Knee Length (AD) Closed Toe Stocking
9291	Medium Compression Thigh Length (AG) Closed Toe Stocking
9280	Medium Compression Knee Length (AD) Open Toe Stocking
9290	Medium Compression Thigh Length (AG) Open Toe Stocking

## VARIMED Cotton For Men

### Ccl II MEDIUM COMPRESSION STOCKINGS (23 – 32 mmHg)

CODE	PRODUCT
8281	Medium Compression Knee Length (AD) Closed Toe Stocking
8291	Medium Compression Thigh Length (AG) Closed Toe Stocking
8280	Medium Compression Knee Length (AD) Open Toe Stocking
8290	Medium Compression Thigh Length (AG) Open Toe Stocking

## VARIMED Classic For Men

### Ccl III STRONG COMPRESSION STOCKINGS (34 -46 mmHg)

CODE	PRODUCT
9380	Strong Compression Knee Length (AD) Open Toe Stocking
9390	Strong Compression Thigh Length (AG) Open Toe Stocking


## VARIMED Cotton For Men

### Ccl III STRONG COMPRESSION STOCKINGS (34 -46 mmHg)

CODE	PRODUCT
8380	Strong Compression Knee Length (AD) Open Toe Stocking
8390	Strong Compression Thigh Length (AG) Open Toe Stocking

## ■ SIZE TABLE

	SIZE						
Leg Circumference (cm)	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cB1	23-27	24-29	26-32	29-35	31-37	33-39	35-41
cC	28-34	30-37	33-40	35-43	37-45	39-47	41-49
cD	27-33	29-36	32-39	34-42	36-44	38-46	40-48
cE	32-38	34-41	37-44	39-47	41-49	43-51	45-53
cF	37-47	40-52	43-56	46-60	50-65	53-68	56-71
cG	43-57	46-62	49-67	52-72	56-76	60-79	64-82



Length of Leg (cm)	
lG Normal	71-82
lG Short	62-70
lD Normal	38-43
lD Short	33-37

**Choosing the Size:** Two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg (lD is for knee-length stockings, lF is for mid-thigh length stockings, lG is for thigh-length and pantyhose stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.



*Knee Length Stocking*

VARIMED Strong



# ■ VARİMED Strong MEDICAL COMPRESSION STOCKING

**Strongest and robust round knit compression stocking**

## **Features:**

- CCL 1,2,3
- Maximum compression stability and stiffness
- High working compression for a deep functional effect
- Perfect fit thanks to soft stretch material, manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Opaque appearance
- Skin friendly - Tested to Oeko-Tex Standard 100

## **Indications:**

- Varicosis in pregnancy with tendency for oedema
- Inflammatory venous disease with tendency for moderate to severe oedema
- Varicosis with tendency for oedema
- Chronic venous insufficiency with tendency for oedema
- Deep venous thrombosis, post-thrombotic syndrome where there is a tendency to oedema
- Lipoedema, stage 1
- Lymphoedema, stage 1



LYCRA® elyaf INVISTA'nın ticari markasıdır.  
LYCRA® fibre is a registered trademark of INVISTA.

## ■ COLORS



**Beige**







*Knee Length Stocking*



*Thigh Length Stocking With Waist Belt*



*Thigh Length Stocking*

## ■ PRODUCT TYPES AND SIZES FOR VARIMED Strong

	STOCKING TYPE					FOOT		LENGTH			SIZE						
Compression Class	AD	AG	AG/T	AT	ATU	OPEN	CLOSED	SHORT	NORMAL	LONG	1	2	3	4	5	6	7
I (18 - 21 mmHg)	●	●	●	●	○	●	●	●	●	●	●	●	●	●	●	●	●
II (23 - 32 mmHg)	●	●	●	●	○	●	●	●	●	●	●	●	●	●	●	●	●
III (34 - 46 mmHg)	●	●	●	●	○	●	○	●	●	●	●	●	●	●	●	●	●

**AD:** KNEE LENGTH COMPRESSION STOCKINGS

**AG:** THIGH LENGTH COMPRESSION STOCKINGS

**AG/T:** THIGH LENGTH WITH WAIST BELT COMPRESSION STOCKINGS

**AT:** PANTY HOSE COMPRESSION STOCKINGS

**ATU:** MATERNITY PANTY HOSE COMPRESSION STOCKINGS

● Available

○ Not Available

## ■ CODES AND DESCRIPTIONS OF VARIMED Strong

### Ccl I LIGHT COMPRESSION STOCKINGS (18 – 21 mmHg)

CODE	PRODUCT
<b>5111</b>	Light Compression Knee Length (AD) Closed Toe Stocking
<b>5131</b>	Light Compression Thigh Length (AG) Closed Toe Stocking
<b>51R1</b>	Light Compression Thigh Length with Waist Belt Right Leg (AG/T) Closed Toe Stocking
<b>51L1</b>	Light Compression Thigh Length with Waist Belt Left Leg (AG/T) Closed Toe Stocking
<b>5141</b>	Light Compression Panty Hose (AT) Closed Toe Stocking
<b>5110</b>	Light Compression Knee Length (AD) Open Toe Stocking
<b>5130</b>	Light Compression Thigh Length (AG) Open Toe Stocking
<b>51R0</b>	Light Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
<b>51L0</b>	Light Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
<b>5140</b>	Light Compression Panty Hose (AT) Open Toe Stocking

## CcI II MEDIUM COMPRESSION STOCKINGS (23 – 32 mmHg)


CODE	PRODUCT
5211	Medium Compression Knee Length (AD) Closed Toe Stocking
5231	Medium Compression Thigh Length (AG) Closed Toe Stocking
52R1	Medium Compression Thigh Length with Waist Belt Right Leg (AG/T) Closed Toe Stocking
52L1	Medium Compression Thigh Length with Waist Belt Left Leg (AG/T) Closed Toe Stocking
5241	Medium Compression Panty Hose (AT) Closed Toe Stocking
5210	Medium Compression Knee Length (AD) Open Toe Stocking
5230	Medium Compression Thigh Length (AG) Open Toe Stocking
52R0	Medium Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
52L0	Medium Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
5240	Medium Compression Panty Hose (AT) Open Toe Stocking

## CcI III STRONG COMPRESSION STOCKINGS (34 -46 mmHg)

CODE	PRODUCT
5310	Strong Compression Knee Length (AD) Open Toe Stocking
5330	Strong Compression Thigh Length (AG) Open Toe Stocking
53R0	Strong Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
53L0	Strong Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
5340	Strong Compression Panty Hose (AT) Open Toe Stocking

## ■ SIZE TABLE

	SIZE						
Leg Circumference (cm)	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cB1	23-27	24-29	26-32	29-35	31-37	33-39	35-41
cC	28-34	30-37	33-40	35-43	37-45	39-47	41-49
cD	27-33	29-36	32-39	34-42	36-44	38-46	40-48
cE	32-38	34-41	37-44	39-47	41-49	43-51	45-53
cF	37-47	40-52	43-56	46-60	50-65	53-68	56-71
cG	43-57	46-62	49-67	52-72	56-76	60-79	64-82



Length of Leg (cm)	
lG Normal	71-82
lG Short	62-70
lD Normal	38-43
lD Short	33-37

**Choosing the Size:** Two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg (lD is for knee-length stockings, lF is for mid-thigh length stockings, lG is for thigh-length and pantyhose stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.

## VARIMED LEG ULCER KIT





# ■ VARIMED LEG ULCER KIT

Varimed Leg Ulcer Kit uses an liner and overstocking combination to create a unique pressure profile for the effective treatment of venous leg ulcers.

## Features:

### Liner Stocking:

Varimed Liner Stocking can be worn 24 hours and has also a therapeutic effect during the night. At night wearing the liner stocking helps to reduce pain and discomfort, helping patient to rest.

- 11-14 mmHg
- Optimum fixing for the wound dressing
- Knitted with antimicrobial, antistatic, therapeutic, anti-odour and thermodynamic Xstatic silver yarn
- Can be ordered separately
- Skin friendly - Tested to Oeko-Tex Standard 100

### Overstocking

- 23-32 mmHg
- Over stocking on the ankle reinforce the compression during the day
- Easy to wear with zip design
- Soft and elastic heel
- High wearing comfort and high efficiency
- Skin friendly - Tested to Oeko-Tex Standard 100

### Combined Stocking System

- 34-46 mmHg
- Ideal combination of liner and overstocking for optimum therapy
- Perfect sizing. 7 sizes

## Indications:

- Venous leg ulcer with oedema
- Prevention of recurrence following healing of venous



LYCRA® elyaf INVISTA'nın ticari markasıdır.  
LYCRA® fibre is a registered trademark of INVISTA.

## ■ COLORS



*Beige*

*Silver*



**Overstocking**



**Liner Stocking**




**Combined Stocking**

## PRODUCT TYPES AND SIZES FOR VARIMED LEG ULCER KIT

	STOCKING TYPE	FOOT		LENGTH	SIZE						
Compression Class	AD	OPEN	CLOSED	NORMAL	1	2	3	4	5	6	7
Liner (11 - 14mmHg)	●	○	●	●	●	●	●	●	●	●	●
Overstocking (23 - 32 mmHg)	●	●	○	●	●	●	●	●	●	●	●
Combined Stocking (34 - 46 mmHg)	●	○	●	●	●	●	●	●	●	●	●

## SIZE TABLE

	SIZE						
Leg Circumference (cm)	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cB1	23-27	24-29	26-32	29-35	31-37	33-39	35-41
cC	28-34	30-37	33-40	35-43	37-45	39-47	41-49
cD	27-33	29-36	32-39	34-42	36-44	38-46	40-48



Length of Leg (cm)	
lD Normal	38-43

**Choosing the Size:** Two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg (lD). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.

# VARIMED SUPPORT STOCKINGS





## SUPPORT STOCKINGS TO PREVENT VENOUS DISEASES



Swollen feet or heavy legs are not only women experience but also recognise the symptoms after standing or sitting at work all day or after a long journey by plane or by car.

Long periods of standing or sitting can trigger blood congestion in the veins, leading to a gradual enlargement of the vein walls. Support Stockings with the correct pressure distribution that decreases from the ankle upwards promote the flow of blood and are therefore an ideal preventative measure for both men and women. Support stocking with correct compression distribution decreasing from the foot upwards accelerate blood flow and are therefore an ideal preventative measure for both men and women.

Support stockings are worn by health persons as a preventive measure against venous problems. They are thinner than compression stockings and pressure on the leg is not as high as with compression stockings. Nevertheless, they actively support the blood flow in the legs. Venous insufficiency is a chronic disease triggered by many risk factors such as sedentary lifestyles, inheritance, excessive weight, pregnancy, hormone treatments. Support stockings have a significant benefit in protecting against these risk factors.

While support stockings can be purchased without prescription, medical compression stockings are generally prescribed by a physician for therapeutic purposes in the case of a venous disease.

### Four compression degrees available for Varimed Support Stockings:

- 40 den 6-9 mmHg (extra light compression)
- 70 den 9-12 mmHg (light compression)
- 140 den 15-20 mmHg (medium compression)
- 280 den 20-25 mmHg (strong compression)

The denier indicate the size (diameter) of the yarn used to weave the stockings. As the yarn denier increases, the compression of the fabric increases, but the socks are less transparent.



VARIMED Glamour 40 Den



## ■ VARIMED Glamour 40 Den

6-9 mmHg

### Features:

- Transparent look
- Perfect fitting
- Soft and smooth touch
- Elegant appearance
- Attractive lace fastening band
- Anatomic and reinforced heel
- Flexible soft toe
- Decreasing support compression from the ankle upwards
- Skin friendly - Tested to Oeko-Tex Standard 100

### Recommended for:

- Ideal for those who stand many hours or sedentary life
- Prevention of minor symptoms of venous insufficiency such as heavy legs and fatigue
- Spider veins



## ■ COLORS



*Skin*

*Bronze*



*Navy*

*Black*



**Panty Hose Stocking**

VARIMED Glamour 70 Den



## ■ VARIMED Glamour 70 Den

9-12 mmHg

### Features:

- Transparent look
- Perfect fitting
- Soft and smooth touch
- Elegant appearance
- Attractive lace fastening band
- Anatomic and reinforced heel
- Flexible soft toe
- Decreasing support compression from the ankle upwards
- Produced by skin friendly – Tested to Oeko -Tex Standard 100

### Recommended for:

- Ideal for those who stand many hours or sedentary life
- Prevention of mild symptoms of venous insufficiency in presence of swollen ankles and feet
- Excites blood circulation and reduce fatigue of your legs
- Spider veins



## ■ COLORS



**Skin**



**Bronze**



**Black**



**Knee Length Stocking**

**VARIMED Glamour 140 Den**





## ■ VARIMED Glamour 140 Den

15-20 mmHg

### Features:

- Slight Transparent look
- Perfect fitting
- Soft and smooth touch
- Elegant appearance
- Attractive lace fastening band
- Anatomic and reinforced heel
- Flexible soft toe
- Decreasing support compression from the ankle upwards
- Skin friendly - Tested to Oeko-Tex Standard 100

### Recommended for:

- Ideal for those who stand many hours or sedentary life
- Prevention of mild symptoms of venous insufficiency in presence of swollen ankles and feet
- Excites blood circulation and reduce fatigue of your legs
- Spider veins
- Reticular varicose veins



## ■ COLORS



*Skin*



*Bronze*



*Black*



**Thigh Length Stocking**

## ■ PRODUCT TYPES AND SIZES FOR VARIMED Glamour

	STOCKING TYPE				FOOT		LENGTH	SIZE					
Compression Class	AD	AG	AT	ATU	OPEN	CLOSED	NORMAL	XS	S	M	L	XL	XXL
40 DEN (6-9 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●
70 DEN (9-12 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●
140 DEN (15-20 mmHg)	●	●	●	○	●	●	●	●	●	●	●	●	●

## ■ CODES AND DESCRIPTIONS OF VARIMED Glamour

### 40 DEN SUPPORT STOCKINGS (6-9 mmHg)

CODE	PRODUCT
2011	40 DEN Glamour Knee Length (AD) Closed Toe Stocking
2031	40 DEN Glamour Thigh Length (AG) Closed Toe Stocking
2041	40 DEN Glamour Panty Hose (AT) Closed Toe Stocking

### 70 DEN SUPPORT STOCKINGS (9-12 mmHg)

CODE	PRODUCT
2111	70 DEN Glamour Knee Length (AD) Closed Toe Stocking
2131	70 DEN Glamour Thigh Length (AG) Closed Toe Stocking
2141	70 DEN Glamour Panty Hose (AT) Closed Toe Stocking

### 140 DEN SUPPORT STOCKINGS (15-20 mmHg)

CODE	PRODUCT
2211	140 DEN Glamour Knee Length (AD) Closed Toe Stocking
2231	140 DEN Glamour Thigh Length (AG) Closed Toe Stocking
2241	140 DEN Glamour Panty Hose (AT) Closed Toe Stocking

## ■ SIZE TABLE FOR VARIMED Glamour 40 DEN,70 DEN and 140 DEN

Below Knee						
AD	XS	S	M	L	XL	XXL
Shoe Size	35-37	36-38	38-39	39-40	40-41	41->
Ankle Circumference (cm)	17-19	18-20	20-23	23-26	26-29	29-32

Thigh Length						
AG	XS	S	M	L	XL	XXL
Thigh Circumference(cm)	43-49	45-51	48-56	50-60	52-62	55-65
Weight in kg	40-48	45-51	52-60	61-69	70-78	79-90
Ankle Circumference(cm)	17-19	18-20	20-23	23-26	26-29	29-32

Panty Hose						
AT	XS	S	M	L	XL	XXL
Hips Circumference(cm)	80-90	85-95	90-100	90-115	95-125	100-130
Weight in kg	40-48	45-51	52-60	61-69	70-78	79-90
Ankle Circumference(cm)	17-19	18-20	20-23	23-26	26-29	29-32



**Knee Length Stocking**



**Panty Hose Stocking**

VARIMED Style 280 Den



## ■ VARIMED Style 280 Den

20-25 mmHg

### Features:

- Support complies with regulations for medical compression stockings
- Strong support
- Perfect fitting
- Soft and smooth touch
- Elegant appearance
- Anatomic and reinforced heel
- Flexible soft toe
- Available with open toe
- Decreasing support compression from the ankle upwards
- Skin friendly – Tested to Oeko-Tex Standard 100

### Recommended for:

- Prevention of oedema from long periods of standing or setendary life,
- Prevention of travel DVT,
- Symptomatic treatment of the initial stages of venous insufficiency
- Excites blood circulation and reduce fatigue of your legs
- Spider veins
- Reticular varicose veins



## ■ COLORS



**Skin**



**Bronze**



**Black**



**Panty Hose Stocking**





***Knee Length Stocking***



***Open Toe***



***Thigh Length Stocking***



***Flat Silicone Topband***

## ■ PRODUCT TYPES AND SIZES FOR VARIMED Style

Compression Class	STOCKING TYPE				FOOT		LENGTH	SIZE						
	AD	AG	AT	ATU	OPEN	CLOSED	NORMAL	1	2	3	4	5	6	7
280 DEN (20-25 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●


## ■ CODES AND DESCRIPTIONS OF VARIMED Style

### 280 DEN SUPPORT STOCKINGS (20-25 mmHg)

CODE	PRODUCT
<b>2311</b>	280 DEN Style Knee Length (AD) Closed Toe Stocking
<b>2331</b>	280 DEN Style Thigh Length (AG) Closed Toe Stocking
<b>2341</b>	280 DEN Style Panty Hose (AT) Closed Toe Stocking
<b>2310</b>	280 DEN Style Knee Length (AD) Open Toe Stocking
<b>2330</b>	280 DEN Style Thigh Length (AG) Open Toe Stocking
<b>2340</b>	280 DEN Style Panty Hose (AT) Open Toe Stocking

## ■ SIZE TABLE

Leg Circumference (cm)	SIZE						
	1	2	3	4	5	6	7
<b>cB</b>	18-20	20-22	22-24	24-26	26-28	28-30	30-32
<b>cB1</b>	23-27	24-29	26-32	29-35	31-37	33-39	35-41
<b>cC</b>	28-34	30-37	33-40	35-43	37-45	39-47	41-49
<b>cD</b>	27-33	29-36	32-39	34-42	36-44	38-46	40-48
<b>cE</b>	32-38	34-41	37-44	39-47	41-49	43-51	45-53
<b>cF</b>	37-47	40-52	43-56	46-60	50-65	53-68	56-71
<b>cG</b>	43-57	46-62	49-67	52-72	56-76	60-79	64-82



**Choosing the Size:** Two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg (lD is for knee-length stockings, lF is for mid-thigh length stockings, lG is for thigh-length and pantyhose stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.

# VARIMED TRAVEL STOCKING



## ■ THROMBOSIS IN THE TRAVEL

Thrombosis can even occur in healthy people who sit for long periods, for instance, during lengthy drives or plane flights.

Wearing travel stockings (support hose) continuously makes the veins narrow slightly and the blood flow more rapidly, making clotting less likely.

Clinical studies have shown that;

Varimed travel stockings are able to reduce the incidence of asymptomatic deep vein thrombosis and of leg oedema. Very large number of passengers would be needed to assess the outcome concerning symptomatic DVT, pulmonary embolism and death.

## ■ VARIMED TRAVEL STOCKING

Varimed travel stockings are designed for the health and comfort of your legs in long-term travels. Long wear-life, perfectly wrapping the legs of stocking, easily wearing and removing the stockings are enabled by double covered yarn which has excellent elasticity. Using all the materials that have Oeko-Tex 100 certificate in production minimizes allergic effects on the skin during wearing of the stocking. Varimed travelling stockings have graduated profile with maximum compression at the ankle, decreasing towards the knee.

### Features:

- High wearing comfort and efficacy
- Soft touch
- Opaque appearance
- Extremely durable yarns
- Easy put on and take off
- Freedom of movements
- Soft, elastic toe and heel
- Perfect fit
- Five sizes
- Skin friendly - Tested to Oeko-Tex Standard 100
- Unisex



### Indications:

In long-term travels;

- Protects against the risk of deep vein thrombosis,
- Accelerates blood flow in the legs by providing graduated compression,
- Prevents swelling of the legs,
- Provides a healthy and comfortable journey.

### Contraindications:

- Decompensated heart failure,
- Advanced peripheral arterial occlusive disease,
- Septic phlebitis,
- Sensory disorders of the limb,
- Advanced peripheral neuropathy,
- Primary chronic arthritis.


## PRODUCT CODES AND DESCRIPTION FOR VARIMED Travel Stockings

CODE	PRODUCT
9161	Travel Stocking
91S1	Antibacterial Travel Stocking with Silver Yarn

## SIZE TABLE

SIZE					
Circumference (cm)	XS	S	M	L	XL
cB	18-20	20-23	23-26	26-29	29-32

**Choosing the Size:** For Varimed Travel Stockings, it is sufficient to know the ankle circumference (cB) measurement to choose the exact size. According to the measurements given in the table, select the appropriate stocking.



## COLORS



*Black*







# VARIMED ANTI-EMBOLISM STOCKINGS



# ■ VARIMED ANTI-EMBOLISM STOCKINGS



## COMPRESSION TREATMENT IN PREVENTING DEEP VEIN THROMBOSIS

Deep vein thrombosis is blood clotting in the deep veins.

A clot that forms in a blood vessel is called a thrombus. Although thrombi can occur in either the superficial or deep leg veins, only those in the deep veins are potentially dangerous. Deep vein thrombosis is dangerous because all or part of the thrombus can break loose, float along in the bloodstream and lodge in a narrow artery in the lung, obstructing blood flow. A moving thrombus is called an embolus.

Deep vein thrombosis shouldn't be confused with phlebitis in varicose veins which is painful but comparatively harmless.

Three factors can contribute to deep vein thrombosis: injury to the lining of the vein; an increased tendency for blood to clot, as can happen with some cancers and rarely with oral contraceptive use; and slowing of the blood flow in the veins, as happens during prolonged bed rest because the calf muscles aren't contracting and squeezing the blood toward the heart. For example, deep vein thrombosis can occur in the heart attack patients who lie in hospital beds for several days with little leg movement or in paraplegics who sit for long periods and whose muscles don't function. Injury or major surgery also can increase the tendency for blood to clot. Thrombosis can even occur in healthy people who sit for long periods, for instance, during lengthy drives or plane flights.

About half the people with deep vein thrombosis have no symptoms at all. In these people, chest pain caused by pulmonary embolism may be the first indication that something is wrong. When deep vein thrombosis causes substantial inflammation and blood flow obstruction, the calf swells and may be painful, tender to the touch and warm. The ankle, foot or thigh may also swell, depending on which veins are involved.

Some thrombi heal by being converted to scar tissue which may damage the valves in the veins. The resulting accumulation of fluid (oedema) can make the ankle swell. The oedema can extend up the leg and even affect the thigh, if the obstruction is high enough in the vein. Oedema is worse toward the end of the day because of the effect of gravity when standing or sitting. Overnight the oedema subsides because the veins empty well when the legs are horizontal.

Although the risk of deep vein thrombosis can't be entirely eliminated, it can be reduced in several ways. One of these ways is wearing anti-embolism stockings continuously makes the veins narrow slightly and the blood flow more rapidly, making clotting less likely. They exert an external pressure which is greater at ankle and reduces at the calf and thigh, this increasing blood velocity within the deep venous system.

### Indications:

- Prophylaxis against deep venous thrombosis in medical and surgical patients.
- Improve resolution of varicose veins by providing graduated compression and increased venous blood flow return.
- Improve blood flow of other venous disorders, such as chronic venous stasis.

### Contraindications:

- Leg conditions that may be exacerbated by stockings include;
  - a) Cellulitis,
  - b) Dermatitis,
  - c) Gangrene,
  - d) Fasciitis,
  - e) Recent skin grafting,
  - f) Panniculitis,
- Severe arteriosclerosis or other ischemic vascular disease.
- Massive oedema of legs or pulmonary oedema from congestive heart failure.
- Extreme deformity of leg



***Thigh Length Stocking***



***Knee Length Stocking***



***Anti-Embolism Stocking with Upper Inspection Hole***



***Anti-Embolism Stocking with Open Toe***

## ■ VARIMED Classic ANTI-EMBOLISM STOCKINGS

### Features:

Varimed Classic Line Anti Embolism Stockings are designed in order to obtain maximum comfort and benefit to patients. Anti-embolism stocking which are produced as knee-length, thigh length and thigh length with waist belt types, has a measurement system that covers wide measurement range. Also, different colours are given to the heel region of each size of stockings to provide easy separation. Stocking has an inspection hole which controls the blood circulation under the foot fingers. Thigh length anti-embolism stocking has dot silicone topband which prevents the stocking from slipping and not distorts the compression profile. Also, an elastic tape which doesn't chafe the leg with excellent elasticity and softness is used during the sewing of thigh length with waist belt anti-embolism stocking.

- 18 mmHg, 23 mmHg, 35 mmHg
- White colour
- Colour coded heel
- High wearing comfort and efficacy
- Soft touch
- Extremely durable yarns
- Easy put on and take off
- Freedom of movements
- Soft, elastic toe and heel
- Perfect sizing. Six sizes and two length
- Skin friendly - Tested to Oeko-Tex Standard 100
- With upper inspection hole and open toe are produced upon request



## ■ PRODUCT TYPES AND SIZES FOR VARIMED Classic ANTI-EMBOLISM STOCKINGS

Compression Class	STOCKING TYPE			FOOT			LENGTH			SIZE					
	AD	AG	AG/T	UPPER INSPECTION HOLE	LOWER INSPECTION HOLE	OPEN	SHORT	NORMAL	LONG	XS	S	M	L	XL	XXL
I (18 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
II (23 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
III (35 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●



# CODES AND DESCRIPTIONS OF VARİMED Classic ANTI-EMBOLISM STOCKINGS

## Ccl I LIGHT COMPRESSION VARİMED CLASSIC ANTI-EMBOLISM STOCKINGS (18 mmHg)

CODE	PRODUCT
<b>6110</b>	Light Compression Knee Length (AD) Upper Inspection Hole Anti-Embolism Stocking
<b>6120</b>	Light Compression Thigh Length (AG) Upper Inspection Hole Anti-Embolism Stocking
<b>6130</b>	Light Compression Thigh Length with Waist Belt (AG/T) Upper Inspection Hole Anti-Embolism Stocking
<b>6111</b>	Light Compression Knee Length (AD) Lower Inspection Hole Anti-Embolism Stocking
<b>6121</b>	Light Compression Thigh Length (AG) Lower Inspection Hole Anti-Embolism Stocking
<b>6131</b>	Light Compression Thigh Length with Waist Belt (AG/T) Lower Inspection Hole Anti-Embolism Stocking
<b>6112</b>	Light Compression Knee Length (AD) Open Toe Anti-Embolism Stocking
<b>6122</b>	Light Compression Thigh Length (AG) Open Toe Anti-Embolism Stocking
<b>6132</b>	Light Compression Thigh Length with Waist Belt (AG/T) Open Toe Anti-Embolism Stocking

## Ccl II MEDIUM COMPRESSION VARİMED CLASSIC ANTI-EMBOLISM STOCKINGS (23 mmHg)


CODE	PRODUCT
<b>6210</b>	Medium Compression Knee Length (AD) Upper Inspection Hole Anti-Embolism Stocking
<b>6220</b>	Medium Compression Thigh Length (AG) Upper Inspection Hole Anti-Embolism Stocking
<b>6230</b>	Medium Compression Thigh Length with Waist Belt (AG/T) Upper Inspection Hole Anti-Embolism Stocking
<b>6211</b>	Medium Compression Knee Length (AD) Lower Inspection Hole Anti-Embolism Stocking
<b>6221</b>	Medium Compression Thigh Length (AG) Lower Inspection Hole Anti-Embolism Stocking
<b>6231</b>	Medium Compression Thigh Length with Waist Belt (AG/T) Lower Inspection Hole Anti-Embolism Stocking
<b>6212</b>	Medium Compression Knee Length (AD) Open Toe Anti-Embolism Stocking
<b>6222</b>	Medium Compression Thigh Length (AG) Open Toe Anti-Embolism Stocking
<b>6232</b>	Medium Compression Thigh Length with Waist Belt (AG/T) Open Toe Anti-Embolism Stocking

## Ccl III STRONG COMPRESSION VARİMED CLASSIC ANTI-EMBOLISM STOCKINGS (35 mmHg)

CODE	PRODUCT
<b>6310</b>	Strong Compression Knee Length (AD) Upper Inspection Hole Anti-Embolism Stocking
<b>6320</b>	Strong Compression Thigh Length (AG) Upper Inspection Hole Anti-Embolism Stocking
<b>6330</b>	Strong Compression Thigh Length with Waist Belt (AG/T) Upper Inspection Hole Anti-Embolism Stocking
<b>6311</b>	Strong Compression Knee Length (AD) Lower Inspection Hole Anti-Embolism Stocking
<b>6321</b>	Strong Compression Thigh Length (AG) Lower Inspection Hole Anti-Embolism Stocking
<b>6331</b>	Strong Compression Thigh Length with Waist Belt (AG/T) Lower Inspection Hole Anti-Embolism Stocking
<b>6312</b>	Strong Compression Knee Length (AD) Open Toe Anti-Embolism Stocking
<b>6322</b>	Strong Compression Thigh Length (AG) Open Toe Anti-Embolism Stocking
<b>6332</b>	Strong Compression Thigh Length with Waist Belt (AG/T) Open Toe Anti-Embolism Stocking

## ■ SIZE TABLE

SIZE						
Circumference (cm)	XS	S	M	L	XL	XXL
cB	18-20	20-23	23-26	26-29	29-32	32-35
cB1	22-27	24-30	27-34	30-37	34-40	38-44
cC	28-34	31-37	34-41	37-44	41-48	44-53
cD	27-33	30-36	33-40	36-43	40-47	43-52
cE	32-38	35-41	38-45	41-48	45-53	49-59
cF	37-47	41-53	47-58	52-62	58-71	64-79
cG	43-57	49-63	53-67	58-73	64-82	72-93
Heel Colour	Beige	Yellow	Green	Red	Blue	Grey



Length of Leg (cm)		
	SHORT	NORMAL
lD	33-37	38-43
lG	62-70	71-82

**Choosing the Size:** For Varimed® Anti-Embolism Stockings two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg ( lD is for knee-length stockings, lG is for thigh-length and thigh-length with waist belt stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.

**Note:** Anti-Embolism Stockings feature the inspection hole under the toes. Models with inspection hole on the upper side and with open toes are produced upon request.

## ■ VARIMED Comfort ANTI-EMBOLISM STOCKINGS

### Features:

Economical solution for patients and hospitals.

Varimed Comfort Line Anti Embolism Stockings are designed in order to obtain maximum comfort and benefit to patients. Anti-embolism stocking which are produced as knee-length, thigh length and thigh length with waist belt types, has a measurement system that covers wide measurement range. Also, different colours are given to the inspection hole of each size of stockings to provide easy separation. Stocking has an inspection hole which controls the blood circulation under the foot fingers. Thigh length anti-embolism stocking has dot silicone topband which prevents the stocking from slipping and not distorts the compression profile. Also, an elastic tape which doesn't chafe the leg with excellent elasticity and softness is used during the sewing of thigh length with waist belt anti-embolism stocking.

- 18 mmHg, 23 mmHg, 35 mmHg
- Ecru colour
- Colour coded inspection hole
- High wearing comfort and efficacy
- Easy put on and take off
- Freedom of movements
- Soft, elastic toe and heel
- Perfect sizing. Six sizes and two length
- Skin friendly - Tested to Oeko-Tex Standard 100
- With upper inspection hole and open toe are produced upon request



## ■ PRODUCT TYPES AND SIZES FOR VARIMED Comfort ANTI-EMBOLISM STOCKINGS

Compression Class	STOCKING TYPE			FOOT			LENGTH			SIZE					
	AD	AG	AG/T	UPPER INSPECTION HOLE	LOWER INSPECTION HOLE	OPEN	SHORT	NORMAL	LONG	XS	S	M	L	XL	XXL
I (18 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
II (23 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
III (35 mmHg)	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●



***Knee Length Stocking***



***Thigh Length Stocking***



***Thigh Length Stockings With Waist Belt***

# CODES AND DESCRIPTIONS OF VARIMED Comfort ANTI-EMBOLISM STOCKINGS

## Ccl I LIGHT COMPRESSION VARIMED COMFORT ANTI-EMBOLISM STOCKINGS (18 mmHg)

CODE	PRODUCT
<b>6110E</b>	Light Compression Knee Length (AD) Upper Inspection Hole Anti-Embolism Stocking
<b>6120E</b>	Light Compression Thigh Length (AG) Upper Inspection Hole Anti-Embolism Stocking
<b>6130E</b>	Light Compression Thigh Length with Waist Belt (AG/T) Upper Inspection Hole Anti-Embolism Stocking
<b>6111E</b>	Light Compression Knee Length (AD) Lower Inspection Hole Anti-Embolism Stocking
<b>6121E</b>	Light Compression Thigh Length (AG) Lower Inspection Hole Anti-Embolism Stocking
<b>6131E</b>	Light Compression Thigh Length with Waist Belt (AG/T) Lower Inspection Hole Anti-Embolism Stocking
<b>6112E</b>	Light Compression Knee Length (AD) Open Toe Anti-Embolism Stocking
<b>6122E</b>	Light Compression Thigh Length (AG) Open Toe Anti-Embolism Stocking
<b>6132E</b>	Light Compression Thigh Length with Waist Belt (AG/T) Open Toe Anti-Embolism Stocking

## Ccl II MEDIUM COMPRESSION VARIMED COMFORT ANTI-EMBOLISM STOCKINGS (23 mmHg)

CODE	PRODUCT
<b>6210E</b>	Medium Compression Knee Length (AD) Upper Inspection Hole Anti-Embolism Stocking
<b>6220E</b>	Medium Compression Thigh Length (AG) Upper Inspection Hole Anti-Embolism Stocking
<b>6230E</b>	Medium Compression Thigh Length with Waist Belt (AG/T) Upper Inspection Hole Anti-Embolism Stocking
<b>6211E</b>	Medium Compression Knee Length (AD) Lower Inspection Hole Anti-Embolism Stocking
<b>6221E</b>	Medium Compression Thigh Length (AG) Lower Inspection Hole Anti-Embolism Stocking
<b>6231E</b>	Medium Compression Thigh Length with Waist Belt (AG/T) Lower Inspection Hole Anti-Embolism Stocking
<b>6212E</b>	Medium Compression Knee Length (AD) Open Toe Anti-Embolism Stocking
<b>6222E</b>	Medium Compression Thigh Length (AG) Open Toe Anti-Embolism Stocking
<b>6232E</b>	Medium Compression Thigh Length with Waist Belt (AG/T) Open Toe Anti-Embolism Stocking


## Ccl III STRONG COMPRESSION VARIMED COMFORT ANTI-EMBOLISM STOCKINGS (35 mmHg)

CODE	PRODUCT
<b>6310E</b>	Strong Compression Knee Length (AD) Upper Inspection Hole Anti-Embolism Stocking
<b>6320E</b>	Strong Compression Thigh Length (AG) Upper Inspection Hole Anti-Embolism Stocking
<b>6330E</b>	Strong Compression Thigh Length with Waist Belt (AG/T) Upper Inspection Hole Anti-Embolism Stocking
<b>6311E</b>	Strong Compression Knee Length (AD) Lower Inspection Hole Anti-Embolism Stocking
<b>6321E</b>	Strong Compression Thigh Length (AG) Lower Inspection Hole Anti-Embolism Stocking
<b>6331E</b>	Strong Compression Thigh Length with Waist Belt (AG/T) Lower Inspection Hole Anti-Embolism Stocking
<b>6312E</b>	Strong Compression Knee Length (AD) Open Toe Anti-Embolism Stocking
<b>6322E</b>	Strong Compression Thigh Length (AG) Open Toe Anti-Embolism Stocking
<b>6332E</b>	Strong Compression Thigh Length with Waist Belt (AG/T) Open Toe Anti-Embolism Stocking



## ■ SIZE TABLE

SIZE						
Circumference (cm)	XS	S	M	L	XL	XXL
cB	18-20	20-23	23-26	26-29	29-32	32-35
cB1	22-27	24-30	27-34	30-37	34-40	38-44
cC	28-34	31-37	34-41	37-44	41-48	44-53
cD	27-33	30-36	33-40	36-43	40-47	43-52
cE	32-38	35-41	38-45	41-48	45-53	49-59
cF	37-47	41-53	47-58	52-62	58-71	64-79
cG	43-57	49-63	53-67	58-73	64-82	72-93
Heel Colour	Beige	Yellow	Green	Red	Blue	Grey



Length of Leg (cm)		
	SHORT	NORMAL
lD	33-37	38-43
lG	62-70	71-82

**Choosing the Size:** Anti-Embolism Stockings two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg (lD is for knee-length stockings, lG is for thigh-length and thigh-length with waist belt stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.

**Note:** Anti-Embolism Stockings feature the inspection hole under the toes. Models with inspection hole on the upper side and with open toes are produced upon request.

# VARIMED DIABETIC SOCKS



## VARIMED DIABETIC SOCKS



Special Sock for diabetics with antimicrobial Xstatic silver yarn

- Without compression and with no upper elastic area
- Completely flat tip seam prevents irritation of the toes
- Knitted with cotton and antimicrobial, antistatic, therapeutic, anti-odour and thermodynamic Xstatic silver yarn
- Socket type is available
- Black, Silver, Navy, Grey and Anthracite colors are available
- Skin friendly - Tested to Oeko-Tex Standard 100
- Unisex



**X-static**  
THE SILVER FIBER®



## COLORS

### Diabetic Sock



**Black**

**Navy**

**Grey**

**Anthracite**

### Diabetic Ankle Sock



**Black**

**Silver**



***Diabetic Sock***



***Diabetic Sock***



***Diabetic Ankle Sock***

## ■ PRODUCT TYPES AND SIZES FOR VARIMED DIABETIC SOCKS

	STOCKING TYPE		FOOT	SIZE		
Color	Diabetic Sock	Diabetic Ankle Sock	CLOSED	35/38	39/42	43/45
Black	●	●	●	●	●	●
Silver	○	●	●	●	●	●
Navy	●	○	●	●	●	●
Grey	●	○	●	●	●	●
Anthracite	●	○	●	●	●	●

## ■ CODES AND DESCRIPTIONS OF VARIMED DIABETIC SOCKS

CODE	PRODUCT
90D12101	Anti-Bacterial Diabetic Sock Black Size:35/38
90D12201	Anti-Bacterial Diabetic Sock Black Size:39/42
90D12301	Anti-Bacterial Diabetic Sock Black Size:43/45
90D12102	Anti-Bacterial Diabetic Sock Navy Size:35/38
90D12202	Anti-Bacterial Diabetic Sock Navy Size:39/42
90D12302	Anti-Bacterial Diabetic Sock Navy Size:43/45
90D12103	Anti-Bacterial Diabetic Sock Grey Size:35/38
90D12203	Anti-Bacterial Diabetic Sock Grey Size:39/42
90D12303	Anti-Bacterial Diabetic Sock Grey Size:43/45
90D12104	Anti-Bacterial Diabetic Sock Anthracite Size:35/38
90D12204	Anti-Bacterial Diabetic Sock Anthracite Size:39/42
90D12304	Anti-Bacterial Diabetic Sock Anthracite Size:43/45
90E12101	Anti-Bacterial Diabetic Ankle Sock Black Size:35/38
90E12201	Anti-Bacterial Diabetic Ankle Sock Black Size:39/42
90E12301	Anti-Bacterial Diabetic Ankle Sock Black Size:43/45
90E12106	Anti-Bacterial Diabetic Ankle Sock Silver Size:35/38
90E12206	Anti-Bacterial Diabetic Ankle Sock Silver Size:39/42
90E12306	Anti-Bacterial Diabetic Ankle Sock Silver Size:43/45







## LYMPHOLOGY

## VARIMED COMPRESSION TREATMENT IN LYMPHOEDEMA



Lymphedema is swelling caused by interference with the normal drainage of lymph back in to the blood. Rarely, lymphedema is obvious at birth. More often, it appears later in life from either congenital (primary lymphoedema) or acquired causes (secondary lymphoedema).

Congenital lymphedema results from having so few lymphatic vessels that they can't handle all the lymph. The problem almost always affects the legs; rarely, it affects the arms. Women are much more likely than men to have congenital lymphedema.

Acquired lymphedema is more common than congenital lymphedema. It typically appears after major surgical treatment, especially after cancer treatment in which lymph nodes and lymphatic vessels are removed or irritated with x-rays. For example, the arm may become prone to swelling after removal of a cancerous breast and associated lymph nodes.

Scarring of repeatedly infected lymphatic vessels also may cause lymphedema, but this is very uncommon except in infection by the tropical parasite *Filaria*.

In acquired lymphedema, the skin looks healthy but it is puffy or swollen. Pressing the area with a finger doesn't leave an indentation, as it does when swelling caused by an accumulation of fluid (edema) results from inadequate blood flow in the veins. In rare instances, the swollen limb is extremely large and the skin is so thick and ridged that it looks almost like elephant skin (elephantiasis).

Lymphedema has no cure. For people with mild lymphedema, medical compression sleeves or stockings can reduce the swelling.

Compression plays an important part in the treatment of lymphedema in the hand, arm, leg or elsewhere. The compressive knitted fabric increase the pressure on the tissue layers in the arm to improve lymphatic drainage from the regional lymph nodes. This prevents renewed swelling (reoedematisation) of the extremity and stabilizes the results of previous treatment.

The air permeable, skin friendly knit with a pressure gradient that steadily decreases from bottom to up ensures physiological lymph flow. The lymphoedema sleeve moulds itself to the limb.

Flat knitted compression stockings or round knitted compression sleeves with high wall stability are used for lymphoedema and lipoedema. They help to keep the oedema constant after decongestion. Varimed offers different qualities with functional characteristics, according to respective indications.





VARIMED STRONG MEDICAL COMPRESSION STOCKINGS





# **VARIMED Strong MEDICAL COMPRESSION STOCKINGS**

Strongest and robust round knit compression stocking

## **Features:**

- CCL 1,2,3
- Maximum compression stability and stiffness
- High working pressure for a deep functional effect
- High wall stability with thickened special knitting system
- Perfect fit thanks to soft stretch material, manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Opaque appearance
- Skin friendly - Tested to Oeko-Tex Standard 100

## **Indications:**

- Varicosis in pregnancy with tendency for oedema
- Inflammatory venous disease with tendency for oedema
- Varicosis with tendency for oedema
- Chronic venous insufficiency with tendency for oedema
- Deep venous thrombosis, post-thrombotic syndrome where there is a tendency to oedema
- Lipoedema, stage 1
- Lymphoedema, stage 1



## **COLORS**



**Beige**



***Knee Length Stocking***



***Thigh Length Stocking With Waist Belt***



***Thigh Length Stocking***

## ■ PRODUCT TYPES AND SIZES FOR VARIMED Strong

	STOCKING TYPE					FOOT		LENGTH			SIZE						
Compression Class	AD	AG	AG/T	AT	ATU	OPEN	CLOSED	SHORT	NORMAL	LONG	1	2	3	4	5	6	7
I (18 - 21 mmHg)	●	●	●	●	○	●	●	●	●	●	●	●	●	●	●	●	●
II (23 - 32 mmHg)	●	●	●	●	○	●	●	●	●	●	●	●	●	●	●	●	●
III (34 - 46 mmHg)	●	●	●	●	○	●	○	●	●	●	●	●	●	●	●	●	●

**AD:** KNEE LENGTH COMPRESSION STOCKINGS

**AG:** THIGH LENGTH COMPRESSION STOCKINGS

**AG/T:** THIGH LENGTH WITH WAIST BELT COMPRESSION STOCKINGS

**AT:** PANTY HOSE COMPRESSION STOCKINGS

**ATU:** MATERNITY PANTY HOSE COMPRESSION STOCKINGS

● Available

○ Not Available

## ■ CODES AND DESCRIPTIONS OF VARIMED Strong

### CcI LIGHT COMPRESSION STOCKINGS (18 – 21 mmHg)

CODE	PRODUCT
<b>5111</b>	Light Compression Knee Length (AD) Closed Toe Stocking
<b>5131</b>	Light Compression Thigh Length (AG) Closed Toe Stocking
<b>51R1</b>	Light Compression Thigh Length with Waist Belt Right Leg (AG/T) Closed Toe Stocking
<b>51L1</b>	Light Compression Thigh Length with Waist Belt Left Leg (AG/T) Closed Toe Stocking
<b>5141</b>	Light Compression Panty Hose (AT) Closed Toe Stocking
<b>5110</b>	Light Compression Knee Length (AD) Open Toe Stocking
<b>5130</b>	Light Compression Thigh Length (AG) Open Toe Stocking
<b>51R0</b>	Light Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
<b>51L0</b>	Light Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
<b>5140</b>	Light Compression Panty Hose (AT) Open Toe Stocking

### **Ccl II MEDIUM COMPRESSION STOCKINGS (23 – 32 mmHg)**


<b>CODE</b>	<b>PRODUCT</b>
<b>5211</b>	Medium Compression Knee Length (AD) Closed Toe Stocking
<b>5231</b>	Medium Compression Thigh Length (AG) Closed Toe Stocking
<b>52R1</b>	Medium Compression Thigh Length with Waist Belt Right Leg (AG/T) Closed Toe Stocking
<b>52L1</b>	Medium Compression Thigh Length with Waist Belt Left Leg (AG/T) Closed Toe Stocking
<b>5241</b>	Medium Compression Panty Hose (AT) Closed Toe Stocking
<b>5210</b>	Medium Compression Knee Length (AD) Open Toe Stocking
<b>5230</b>	Medium Compression Thigh Length (AG) Open Toe Stocking
<b>52R0</b>	Medium Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
<b>52L0</b>	Medium Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
<b>5240</b>	Medium Compression Panty Hose (AT) Open Toe Stocking

### **Ccl III STRONG COMPRESSION STOCKINGS (34 -46 mmHg)**

<b>CODE</b>	<b>PRODUCT</b>
<b>5310</b>	Strong Compression Knee Length (AD) Open Toe Stocking
<b>5330</b>	Strong Compression Thigh Length (AG) Open Toe Stocking
<b>53R0</b>	Strong Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
<b>53L0</b>	Strong Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
<b>5340</b>	Strong Compression Panty Hose (AT) Open Toe Stocking

## ■ SIZE TABLE

	SIZE						
Leg Circumference (cm)	1	2	3	4	5	6	7
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32
cB1	23-27	24-29	26-32	29-35	31-37	33-39	35-41
cC	28-34	30-37	33-40	35-43	37-45	39-47	41-49
cD	27-33	29-36	32-39	34-42	36-44	38-46	40-48
cE	32-38	34-41	37-44	39-47	41-49	43-51	45-53
cF	37-47	40-52	43-56	46-60	50-65	53-68	56-71
cG	43-57	46-62	49-67	52-72	56-76	60-79	64-82



Length of Leg (cm)	
lG Normal	71-82
lG Short	62-70
lD Normal	38-43
lD Short	33-37

**Choosing the Size:** Two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg (lD is for knee-length stockings, lF is for mid-thigh length stockings, lG is for thigh-length and pantyhose stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.



# VARIMED Forte MEDICAL COMPRESSION STOCKING



# VARIMED Forte MEDICAL COMPRESSION STOCKINGS

Flat knit special design for lymphoedema

## Features:

- CCL 1,2,3
- Specially designed for treating lymphatic diseases
- Highly effective working pressure
- High wall stability with thickened special knitting system
- Perfect fit thanks to soft stretch material, manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Anatomic design
- Breathable and moisture regulating
- Robust and durable
- Skin friendly - Tested to Oeko-Tex Standard 100

## Indications:

- Primary and secondary lymphoedema stage 1 and 2
- Lipoedema stages 1 and 2
- Postoperative and posttraumatic oedema
- Severe phlebological conditions with a tendency to develop oedema
- Venus ulcer treatment
- Lymphatic ulcer treatment
- Kleoid treatment, e.g. after burn



## COLORS



**Beige**





***Panty Hose Stocking***



***Knee Length Stocking***



***Thigh Length Stocking With Waist Belt***

## ■ PRODUCT TYPES AND SIZES FOR VARIMED Forte

	STOCKING TYPE					FOOT		LENGTH			SIZE						
Compression Class	AD	AG	AG/T	AT	ATU	OPEN	CLOSED	SHORT	NORMAL	LONG	1	2	3	4	5	6	7
I (18 - 21 mmHg)	●	●	●	●	○	●	●	●	●	●	●	●	●	●	●	●	●
II (23 - 32 mmHg)	●	●	●	●	○	●	●	●	●	●	●	●	●	●	●	●	●
III (34 - 46 mmHg)	●	●	●	●	○	●	○	●	●	●	●	●	●	●	●	●	●

**AD:** KNEE LENGTH COMPRESSION STOCKINGS

**AG:** THIGH LENGTH COMPRESSION STOCKINGS

**AG/T:** THIGH LENGTH WITH WAIST BELT COMPRESSION STOCKINGS

**AT:** PANTY HOSE COMPRESSION STOCKINGS

**ATU:** MATERNITY PANTY HOSE COMPRESSION STOCKINGS

● Available

○ Not Available

## ■ CODES AND DESCRIPTIONS OF VARIMED Forte

### CcI I LIGHT COMPRESSION STOCKINGS (18-21 mmHg)

CODE	PRODUCT
<b>5111D</b>	Light Compression Knee Length (AD) Closed Toe Stocking
<b>5131D</b>	Light Compression Thigh Length (AG) Closed Toe Stocking
<b>51R1D</b>	Light Compression Thigh Length with Waist Belt Right Leg (AG/T) Closed Toe Stocking
<b>51L1D</b>	Light Compression Thigh Length with Waist Belt Left Leg (AG/T) Closed Toe Stocking
<b>5141D</b>	Light Compression Panty Hose (AT) Closed Toe Stocking
<b>5110D</b>	Light Compression Knee Length (AD) Open Toe Stocking
<b>5130D</b>	Light Compression Thigh Length (AG) Open Toe Stocking
<b>51R0D</b>	Light Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
<b>51L0D</b>	Light Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
<b>5140D</b>	Light Compression Panty Hose (AT) Open Toe Stocking

## Ccl II MEDIUM COMPRESSION STOCKINGS (23 – 32 mmHg)

CODE	PRODUCT
<b>5211D</b>	Medium Compression Knee Length (AD) Closed Toe Stocking
<b>5231D</b>	Medium Compression Thigh Length (AG) Closed Toe Stocking
<b>52R1D</b>	Medium Compression Thigh Length with Waist Belt Right Leg (AG/T) Closed Toe Stocking
<b>52L1D</b>	Medium Compression Thigh Length with Waist Belt Left Leg (AG/T) Closed Toe Stocking
<b>5241D</b>	Medium Compression Panty Hose (AT) Closed Toe Stocking
<b>5210D</b>	Medium Compression Knee Length (AD) Open Toe Stocking
<b>5230D</b>	Medium Compression Thigh Length (AG) Open Toe Stocking
<b>52R0D</b>	Medium Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
<b>52L0D</b>	Medium Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
<b>5240D</b>	Medium Compression Panty Hose (AT) Open Toe Stocking

## Ccl III STRONG COMPRESSION STOCKINGS (34 -46 mmHg)


CODE	PRODUCT
<b>5310D</b>	Strong Compression Knee Length (AD) Open Toe Stocking
<b>5330D</b>	Strong Compression Thigh Length (AG) Open Toe Stocking
<b>53R0D</b>	Strong Compression Thigh Length with Waist Belt Right Leg (AG/T) Open Toe Stocking
<b>53L0D</b>	Strong Compression Thigh Length with Waist Belt Left Leg (AG/T) Open Toe Stocking
<b>5340D</b>	Strong Compression Panty Hose (AT) Open Toe Stocking





## ■ SIZE TABLE

	SIZE								
Leg Circumference (cm)	0	1	2	3	4	5	6	7	8
cB	16-18	18-20	20-22	22-24	24-26	26-28	28-30	30-32	32-35
cB1	21-25	23-27	24-29	26-32	29-35	31-37	33-39	35-41	37-44
cC	26-31	28-34	30-37	33-40	35-43	37-45	39-47	41-49	43-52
cD	25-30	27-33	29-36	32-39	34-42	36-44	38-46	40-48	42-51
cE	29-35	32-38	34-41	37-44	39-47	41-49	43-51	45-53	47-55
cF	34-43	37-47	40-52	43-56	46-60	50-65	53-68	56-71	60-78
cG	40-53	43-57	46-62	49-67	52-72	56-76	60-79	64-82	69-92



Length of Leg (cm)	
lG Normal	71-82
lG Short	62-70
lD Normal	38-43
lD Short	33-37

**Choosing the Size:** Two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg (lD is for knee-length stockings, lF is for mid-thigh length stockings, lG is for thigh-length and pantyhose stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.

## VARIMED Forte TOE CAP



## ■ VARIMED Forte TOE CAP

Flat knit special design for lymphoedema

### Features:

- CCL 1,2,3
- Stitching system that does not interfere with lymphatic drainage and provides the comfort of the toes
- High wall stability with thickened special knitting system
- Breathable and moisture regulating
- Robust and durable
- With small toe and without small toe
- 7 sizes
- Skin friendly - Tested to Oeko-Tex Standard 100

### INTEGRATED TOE CAP

- Toe cap is integrated with below knee or thigh length stocking as seamless

### INDIVIDUAL TOE CAP

- Individual care options, e.g. combination with other products

### Indications:

- Lymphoedema of the forefoot and toes
- Lymphoedema of the dorsum of the foot



## ■ COLORS



**Beige**



## ■ CODES AND DESCRIPTIONS OF VARIMED Forte TOE CAP

### CcI I LIGHT COMPRESSION STOCKINGS (18 – 21 mmHg)

CODE	PRODUCT
<b>5102D</b>	Light Compression Toe Cap
<b>5112D</b>	Light Compression Knee Length Stocking with Toe Cap (AD)
<b>5132D</b>	Light Compression Thigh Length Stocking with Toe Cap (AG)

### CcI II MEDIUM COMPRESSION STOCKINGS (23 – 32 mmHg)

CODE	PRODUCT
<b>5202D</b>	Medium Compression Toe Cap Stocking
<b>5212D</b>	Medium Compression Knee Length Stocking with Toe Cap (AD)
<b>5232D</b>	Medium Compression Thigh Length Stocking with Toe Cap (AG)

### CcI III STRONG COMPRESSION STOCKINGS (34 -46 mmHg)

CODE	PRODUCT
<b>5302D</b>	Strong Compression Toe Cap Stocking
<b>5312D</b>	Strong Compression Knee Length Stocking with Toe Cap (AD)
<b>5332D</b>	Strong Compression Thigh Length Stocking with Toe Cap (AG)




***Integrated Toe Cap***



***Individual Toe Cap***

## ■ SIZE TABLE FOR TOE CAP

		SIZE						
Toe Circumference (cm)		1	2	3	4	5	6	7
1.Toe	Top	7,5	8,0	8,5	9,5	10,0	10,5	11,0
	Bottom	7,5	8,0	8,5	9,5	10,0	10,5	11,0
2.Toe	Top	5,0	5,5	6,0	6,5	7,0	7,5	8,0
	Bottom	5,0	5,5	6,0	6,5	7,0	7,5	8,0
3.Toe	Top	4,5	5,0	5,5	6,0	6,5	7,0	7,5
	Bottom	4,5	5,0	5,5	6,0	6,5	7,0	7,5
4.Toe	Top	4,5	5,0	5,5	6,0	6,5	7,0	7,5
	Bottom	4,5	5,0	5,5	6,0	6,5	7,0	7,5
5.Toe	Top	4,5	5,0	5,5	6,0	6,5	7,0	7,5
	Bottom	4,5	5,0	5,5	6,0	6,5	7,0	7,5
cA		18-20	20-22	22-24	24-26	26-28	28-30	30-32



**Choosing the Size:** To choose the exact size of Varimed® Forte Toe Cap, it is sufficient to determine the foot Circumference (cA), The ankle circumference (cB) and length of the leg (lD is for knee-length stockings, lG is for thigh-length). You can choose the suitable size from the measurements in the table. Do not use the toe cap in case of patient's whole leg and finger measurements do not correspond to sizes indicated in table.



# VARIMED Strong LYMPHOEDEMA ARMSLEEVE



# VARIMED Strong LYMPHOEDEMA ARMSLEEVES

Special Design Round Knit For Lymphoedema

## Features:

- CCL 1,2,3
- Specially designed for treating lymphatic diseases
- High wall stability with thickened special knitting system
- Perfect fit thanks to soft stretch material, manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Breathable and moisture regulating
- The upper edge is coated with silicone to ensure adherence to limb
- Support belt version
- Perfect sizing. 7 sizes and 2 length
- Skin friendly - Tested to Oeko-Tex Standard 100

## Indications:

- Primary and secondary stage 1 lymphoedema of the arm
- Shoulder and arm vein thrombosis
- Postoperative and posttraumatic oedema
- Not suitable for pronounced oedema



## COLORS



**Beige**



**Lymphoedema Armsleeve Without Hand**



***Lymphoedema Armsleeve With Hand and Shoulder cap Belt***



***Lymphoedema Armsleeve Without Hand and dot silicone topband***



***Lymphoedema Armsleeve With Hand and dot silicone topband***

## PRODUCT TYPES AND SIZES FOR VARIMED Strong LYMPHOEDEMA ARMSLEEVES

### CODES AND DESCRIPTIONS OF VARIMED Strong LYMPHOEDEMA ARMSLEEVES

#### Ccl I LIGHT COMPRESSION ARMSLEEVES (18 – 21 mmHg)

CODE	PRODUCT
511K	Light Compression Without Hand (CG) Armsleeve
512K	Light Compression With Hand (AG) Armsleeve
513K	Light Compression without Hand with Shouldercap Belt (CG/H) Armsleeve
514K	Light Compression with Hand with Shouldercap Belt (AG/H) Armsleeve
515K	Light Compression Without Hand and Dot Silicone Topband (CG) Armsleeve
516K	Light Compression With Hand and Dot Silicone Topband (AG) Armsleeve

#### Ccl II MEDIUM COMPRESSION ARMSLEEVES (23 – 32 mmHg)

CODE	PRODUCT
521K	Medium Compression Without Hand (CG) Armsleeve
522K	Medium Compression With Hand (AG) Armsleeve
523K	Medium Compression without Hand with Shouldercap Belt (CG/H) Armsleeve
524K	Medium Compression with Hand with Shouldercap Belt (AG/H) Armsleeve
525K	Medium Compression Without Hand and Dot Silicone Topband (CG) Armsleeve
526K	Medium Compression With Hand and Dot Silicone Topband (AG) Armsleeve

#### Ccl III STRONG COMPRESSION ARMSLEEVES (34 -46 mmHg)

CODE	PRODUCT
531K	Strong Compression Without Hand (CG) Armsleeve
532K	Strong Compression With Hand (AG) Armsleeve
533K	Strong Compression without Hand with Shouldercap Belt (CG/H) Armsleeve
534K	Strong Compression with Hand with Shouldercap Belt (AG/H) Armsleeve
535K	Strong Compression Without Hand and Dot Silicone Topband (CG) Armsleeve
536K	Strong Compression With Hand and Dot Silicone Topband (AG) Armsleeve

## VARIMED Strong COTTON LYMPHOEDEMA ARMSLEEVES





## ■ VARIMED Strong COTTON LYMPHOEDEMA ARMSLEEVES

Special design cotton round knit

Varimed Strong Cotton Lymphoedema Armsleeves are comfortably worn by patients with sensitive skin or allergy sufferers due to skin compatibility. The inner surface of the armsleeve that touches the skin is 100% cotton.

### Features:

- CCL1, CCL2, CCL3
- Specially designed for treating lymphatic diseases
- High wall stability with thickened special knitting system
- Perfect fit thanks to soft stretch material, manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Breathable and moisture regulating
- The upper edge is coated with silicone to ensure adherence to limb
- Support belt version
- Perfect sizing. 7 sizes and 2 length
- Skin friendly - Tested to Oeko-Tex Standard 100

### Indications:

- Primary and secondary stage 1 lymphoedema of the hand
- Shoulder and arm vein thrombosis
- Postoperative and posttraumatic oedema
- Not suitable for pronounced oedema



## ■ COLORS



**Beige**



***Lymphoedema Armsleeve Without Hand and dot silicone topband***



***Lymphoedema Armsleeve with Hand and Shouldercap Belt***



***Lymphoedema Armsleeve With Glove and dot silicone topband***

## PRODUCT TYPES AND SIZES FOR VARIMED Strong COTTON LYMPHOEDEMA ARMSLEEVES

### CODES AND DESCRIPTIONS OF VARIMED Strong COTTON LYMPHOEDEMA ARMSLEEVES

#### Ccl I LIGHT COMPRESSION ARMSLEEVES (18 – 21 mmHg)

CODE	PRODUCT
511P	Light Compression Without Hand (CG) Armsleeve
512P	Light Compression With Hand (AG) Armsleeve
512G	Light Compression With Glove (AG)Long Finger Armsleeve
512S	Light Compression With Glove (AG)Short Finger Armsleeve
513P	Light Compression without Hand with Shouldercap Belt (CG/H) Armsleeve
514P	Light Compression with Hand with Shouldercap Belt (AG/H) Armsleeve
514G	Light Compression with Glove with Shouldercap Belt (AG/H)Long Finger Armsleeve
514S	Light Compression with Glove with Shouldercap Belt (AG/H)Short Finger Armsleeve
515P	Light Compression Without Hand and Dot Silicone Topband (CG) Armsleeve
516P	Light Compression With Hand and Dot Silicone Topband (AG) Armsleeve
516G	Light Compression With Glove and Dot Silicone Topband (AG)Long Finger Armsleeve
516S	Light Compression With Glove and Dot Silicone Topband (AG)Short Finger Armsleeve

#### Ccl II MEDIUM COMPRESSION ARMSLEEVES (23 – 32 mmHg)

CODE	PRODUCT
521P	Medium Compression Without Hand (CG) Armsleeve
522P	Medium Compression With Hand (AG) Armsleeve
522G	Medium Compression With Glove( AG)Long Finger Armsleeve
522S	Medium Compression With Glove (AG)Short Finger Armsleeve
523P	Medium Compression without Hand with Shouldercap Belt (CG/H) Armsleeve
524P	Medium Compression with Hand with Shouldercap Belt (AG/H) Armsleeve
524G	Medium Compression with Glove with Shouldercap Belt (AG/H)Long Finger Armsleeve
524S	Medium Compression with Glove with Shouldercap Belt( AG/H)Short Finger Armsleeve
525P	Medium Compression Without Hand and Dot Silicone Topband (CG) Armsleeve
526P	Medium Compression With Hand and Dot Silicone Topband (AG) Armsleeve
526G	Medium Compression With Glove and Dot Silicone Topband (AG)Long Finger Armsleeve
526S	Medium Compression With Glove and Dot Silicone Topband (AG)Short Finger Armsleeve


## Ccl III STRONG COMPRESSION ARMSLEEVES (34 -46 mmHg)

CODE	PRODUCT
<b>531P</b>	Strong Compression Without Hand (CG) Armsleeve
<b>532P</b>	Strong Compression With Hand (AG) Armsleeve
<b>532G</b>	Strong Compression With Glove (AG) Long Finger Armsleeve
<b>532S</b>	Strong Compression With Glove (AG) Short Finger Armsleeve
<b>533P</b>	Strong Compression without Hand with Shouldercap Belt (CG/H) Armsleeve
<b>534P</b>	Strong Compression with Hand with Shouldercap Belt (AG/H) Armsleeve
<b>534G</b>	Strong Compression with Glove with Shouldercap Belt (AG/H)Long Finger Armsleeve
<b>534S</b>	Strong Compression with Glove with Shouldercap Belt (AG/H)Short Finger Armsleeve
<b>535P</b>	Strong Compression Without Hand and Dot Silicone Topband (CG) Armsleeve
<b>536P</b>	Strong Compression With Hand and Dot Silicone Topband (AG) Armsleeve
<b>536G</b>	Strong Compression With Glove and Dot Silicone Topband (AG) Long Finger Armsleeve
<b>536S</b>	Strong Compression With Glove and Dot Silicone Topband (AG) Short Finger Armsleeve

## ■ SIZE TABLE

	SIZE						
Leg Circumference (cm)	1	2	3	4	5	6	7
<b>cA</b>	15-18	16,5-19	18-21	20-23	22-25	24-27	26-29
<b>cC</b>	12-13	13-15	15-17	17-19	19-21	21-23	23-25
<b>cD</b>	16-21	20-24	23-27	26-30	28-33	31-36	34-39
<b>cE</b>	17-22	20-26	23-29	26-32	29-35	32-38	35-41
<b>cG</b>	19-25	22-29	25-32	29-36	33-40	37-44	41-48

Length (cm)	Short	Normal
<b>lG</b>	Up to 44 (short<44)	Over 44 cm (normal>44)



**Choosing the size:** To choose the exact size of Varimed® Lymphoedema Armsleeves it is sufficient to determine the sizes of wrist (cC), elbow (cE), topmost part of the arm (cG) and length of arm (lG). Further measurements are provided to verify the size of the armsleeve. You can choose the suitable size from the measurements in the table. Do not use the armsleeve in case of patient's whole arm measurements do not correspond to sizes indicated in table.





# VARIMED Strong COTTON LYMPHOEDEMA GLOVES



## ■ VARIMED Strong COTTON LYMPHOEDEMA GLOVES

Special design cotton round knit

Varimed Strong Cotton Lymphoedema Armsleeves are comfortably worn by patients with sensitive skin or allergy sufferers due to skin compatibility. The inner surface of the armsleeve that touches the skin is 100% cotton.

### Features:

- CCL1,2,3
- Specially designed for treating lymphatic diseases
- High wall stability with thickened special knitting system
- Perfect fit thanks to soft stretch material, manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Breathable and moisture regulating
- 7 sizes
- Short and long finger type
- Individual care options, e.g. combination with other products
- Skin friendly - Tested to Oeko-Tex Standard 100

### Indications:

- Primary and secondary stage 1 lymphoedema of the hand
- Prevention of accumulation of lymph fluid in the hand area while wearing a medical compression armsleeve
- Postoperative and posttraumatic oedema
- Not suitable for pronounced oedema



## ■ COLORS



**Beige**



***Lymphoedema Glove with Short Finger***



***Lymphoedema Glove with Long Finger***

# PRODUCT TYPES AND SIZES FOR VARIMED Strong COTTON LYMPHOEDEMA GLOVES

## CODES AND DESCRIPTIONS OF THE PRODUCTS

### CcI LIGHT COMPRESSION STOCKINGS (18 – 21 mmHg)

CODE	PRODUCT
517P	Light Compression Glove with Short Finger
518P	Light Compression Glove with Long Finger

### CcI II MEDIUM COMPRESSION GLOVES (23 – 32 mmHg)

CODE	PRODUCT
527P	Medium Compression Glove with Short Finger
528P	Medium Compression Glove with Long Finger

### CcI III STRONG COMPRESSION GLOVES (34 -46 mmHg)

CODE	PRODUCT
537P	Strong Compression Glove with Short Finger
538P	Strong Compression Glove with Long Finger

# ■ SIZE TABLE


## Circumference of Hand(cm)

	SIZE						
Hand Circumference (cm)	1	2	3	4	5	6	7
cA	15-18	16,5-19	18-21	20-23	22-25	24-27	26-29
cB	15,5-19	17-20	18,5-22	20,5-24	22,5-26	24,5-28	26,5-30
cC	12-13	13-15	15-17	17-19	19-21	21-23	23-25

Length of Hand(cm)	
lAB	3-5
lAC	9-12

## Circumference of Fingers(cm)

		SIZE								
Finger Circumference (cm)		1	2	3	4	5	6	7	Length of Fingers	
Thumb	Top	5,5	6,0	6,5	7,0	7,5	8,0	9,0	Thumb	3,5-4,5
	Bottom	6,5	7,0	7,5	8,0	8,5	9,0	10,0		
Index Finger	Top	5,0	5,5	6,0	6,5	7,0	7,5	8,5	Index Finger	4,5-5,5
	Bottom	6,0	6,5	7,0	7,5	8,0	8,5	9,5		
Middle Finger	Top	5,0	5,5	6,0	6,5	7,0	7,5	8,5	Middle Finger	5,5-6,5
	Bottom	6,0	6,5	7,0	7,5	8,0	8,5	9,5		
Ring Finger	Top	4,5	5,0	5,5	6,0	6,5	7,0	8,0	Ring Finger	4,5-5,5
	Bottom	5,5	6,0	6,5	7,0	7,5	8,0	9,0		
Little Finger	Top	4,0	4,5	5,0	5,5	6,0	6,5	7,5	Little Finger	3,5-4,5
	Bottom	5,0	5,5	6,0	6,5	7,0	7,5	8,5		



**Choosing the size:** : To choose the exact size of Varimed® Lymphoedema Gloves, it is sufficient to determine the sizes of wrist (cC),hands(cA) and palms of hands elbow (cB).You can choose the suitable size from the measurements in the table. Do not use the glove in case of patient's whole arm and finger measurements do not correspond to sizes indicated in table.







# VARIMED Forte LYMPHOEDEMA ARMSLEEVE



## ■ VARIMED Forte LYMPHOEDEMA ARMSLEEVES

Flat knit special design for lymphoedema

### Features:

- CCL 1,2,3
- Specially designed for treating lymphatic diseases
- Highly effective working pressure
- High wall stability with thickened special knitting system
- Perfect fit thanks to soft stretch material, manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Anatomic design
- Breathable and moisture regulating
- Robust and durable
- The upper edge is coated with silicone to ensure adherence to limb
- Belt version available
- Perfect sizing. 7 sizes and 2 length
- Skin friendly - Tested to Oeko-Tex Standard 100

### Indications:

- Primary and secondary lymphoedema stage 1 and 2
- Lipoedema stages 1 and 2
- Postoperative and posttraumatic oedema
- Severe phlebological conditions with a tendency to develop oedema
- Kleoid treatment, e.g. after burns



## ■ COLORS



**Beige**



## PRODUCT TYPES AND SIZES FOR VARIMED Forte LYMPHOEDEMA ARMSLEEVES

### CODES AND DESCRIPTIONS OF VARIMED Forte LYMPHOEDEMA ARMSLEEVES

#### CcI I LIGHT COMPRESSION ARMSLEEVES (18 - 21 mmHg)

CODE	PRODUCT
<b>513D</b>	Light Compression without Hand with Shouldercap Belt(CG/H) Armsleeve
<b>514D</b>	Light Compression with Hand with Shouldercap Belt(AG/H) Armsleeve
<b>515D</b>	Light Compression Without Hand and Dot Silicone Topband(CG) Armsleeve
<b>516D</b>	Light Compression With Hand and Dot Silicone Topband(AG) Armsleeve

#### CcI II MEDIUM COMPRESSION ARMSLEEVES (23 - 32 mmHg)

CODE	PRODUCT
<b>523D</b>	Medium Compression without Hand with Shouldercap Belt(CG/H) Armsleeve
<b>524D</b>	Medium Compression with Hand with Shouldercap Belt(AG/H) Armsleeve
<b>525D</b>	Medium Compression Without Hand and Dot Silicone Topband(CG) Armsleeve
<b>526D</b>	Medium Compression With Hand and Dot Silicone Topband(AG) Armsleeve


#### CcI III STRONG COMPRESSION ARMSLEEVES (34 - 46 mmHg)

CODE	PRODUCT
<b>533D</b>	Strong Compression without Hand with Shouldercap Belt(CG/H) Armsleeve
<b>534D</b>	Strong Compression with Hand with Shouldercap Belt(AG/H) Armsleeve
<b>535D</b>	Strong Compression Without Hand and Dot Silicone Topband(CG) Armsleeve
<b>536D</b>	Strong Compression With Hand and Dot Silicone Topband(AG) Armsleeve

■ SIZE TABLE

	SIZE						
Circumference (cm)	1	2	3	4	5	6	7
cA	15-18	16,5-19	18-21	20-23	22-25	24-27	26-29
cC	12-13	13-15	15-17	17-19	19-21	21-23	23-25
cD	16-21	20-24	23-27	26-30	28-33	31-36	34-39
cE	17-22	20-26	23-29	26-32	29-35	32-38	35-41
cG	19-25	22-29	25-32	29-36	33-40	37-44	41-48

Length (cm)	Short	Normal
ℓG	Up to 44 (short<44)	Over 44 cm (normal>44)



**Choosing the size:** To choose the exact size of Varimed® Lymphoedema Armsleeves it is sufficient to determine the sizes of wrist (cC), elbow (cE), topmost part of the arm (cG) and length of arm (ℓG). Further measurements are provided to verify the size of the armsleeve. You can choose the suitable size from the measurements in the table. Do not use the armsleeve in case of patient's whole arm measurements do not correspond to sizes indicated in table.





***Lymphoedema Armsleeve Without Hand and Shoulder cap Belt***



***Lymphoedema Armsleeve With Hand  
and dot silicone topband***



***Lymphoedema Armsleeve Without Hand  
and dot silicone topband***



# VARIMED Forte LYMPHOEDEMA GLOVES





## ■ VARIMED Forte LYMPHOEDEMA GLOVES

Flat knit special design for lymphoedema

### Features:

- CCL1, CCL2, CCL3
- Specially designed for treating lymphatic diseases
- High wall stability with thickened special knitting system
- Perfect fit thanks to soft stretch material, manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Anatomic design
- Breathable and moisture regulating
- 7 sizes
- Short and long finger type
- Skin friendly - Tested to Oeko-Tex Standard 100

### Indications:

- Primary and secondary stage 1 and 2 lymphoedema of the hand
- Prevention of accumulation of lymph fluid in the hand area while wearing a medical compression armsleeve
- Postoperative and posttraumatic oedema



## ■ INDIVIDUAL GLOVE

- Individual care options, e.g. combination with other products

## ■ INTEGRATED GLOVE

- Glove is integrated with armsleeve as seamless

## ■ COLORS



**Beige**



***Integrated Glove***



***Individual Glove***



## PRODUCT TYPES AND SIZES FOR VARIMED Forte LYMPHOEDEMA GLOVES

### CODES AND DESCRIPTIONS OF VARIMED Forte LYMPHOEDEMA GLOVES

#### Ccl I LIGHT COMPRESSION ARMSLEEVES WITH GLOVE (18 – 21 mmHg)

CODE	PRODUCT
<b>514DG</b>	Light Compression with Glove with Shouldercap Belt(AG/H)Long Finger Armsleeve
<b>516DG</b>	Light Compression With Glove and Dot Silicone Topband(AG/H)Long Finger Armsleeve

#### Ccl II MEDIUM COMPRESSION ARMSLEEVES WITH GLOVE (23 – 32 mmHg)

CODE	PRODUCT
<b>524DG</b>	Medium Compression with Glove with Shouldercap Belt(AG/H)Long Finger Armsleeve
<b>526DG</b>	Medium Compression With Glove and Dot Silicone Topband(AG/H)Long Finger Armsleeve

#### Ccl III MEDIUM COMPRESSION ARMSLEEVES WITH GLOVE (34 – 46 mmHg)

CODE	PRODUCT
<b>534DG</b>	Strong Compression with Glove with Shouldercap Belt(AG/H)Long Finger Armsleeve
<b>536DG</b>	Strong Compression With Glove and Dot Silicone Topband(AG/H)Long Finger Armsleeve

#### Ccl I LIGHT COMPRESSION GLOVES (18 – 21 mmHg)

CODE	PRODUCT
<b>517D</b>	Light Compression Glove with Short Finger
<b>518D</b>	Light Compression Glove with Long Finger

#### Ccl II MEDIUM COMPRESSION GLOVES (23 – 32 mmHg)

CODE	PRODUCT
<b>527D</b>	Medium Compression Glove with Short Finger
<b>528D</b>	Medium Compression Glove with Long Finger

#### Ccl III STRONG COMPRESSION GLOVES (34 -46 mmHg)

CODE	PRODUCT
<b>537D</b>	Strong Compression Glove with Short Finger
<b>538D</b>	Strong Compression Glove with Long Finger

VARIMED Forte SEAMLESS LYMPHOEDEMA GLOVE





## ■ VARIMED Forte SEAMLESS LYMPHOEDEMA GLOVES

Flat knit and seamless special design for lymphoedema

### Features:

- CCL 1,2,3
- Specially designed for treating lymphatic diseases
- High wall stability with thickened special knitting system
- Perfect fit thanks to soft stretch material, manufactured using medical type Lycra with soft stretch characteristics and the low hysteresis
- Breathable and moisture regulating
- 7 sizes
- Short and long finger type
- Individual care options, e.g. combination with other products
- Skin friendly - Tested to Oeko-Tex Standard 100

### Indications:

- Primary and secondary stage 1 and 2 lymphoedema of the hand
- Prevention of accumulation of lymph fluid in the hand area while wearing a medical compression armsleeve
- Postoperative and posttraumatic oedema



## ■ COLORS



**Beige**

**PRODUCT TYPES AND SIZES FOR VARIMED Forte  
SEAMLESS LYMPHOEDEMA GLOVES**

**CODES AND DESCRIPTIONS FOR VARIMED Forte  
SEAMLESS LYMPHOEDEMA GLOVES**

**CcI I LIGHT COMPRESSION GLOVES (18 – 21 mmHg)**

CODE	PRODUCT
519D	Light Compression Seamless Glove with Long Finger

**CcI II MEDIUM COMPRESSION GLOVES (23 – 32 mmHg)**

CODE	PRODUCT
529D	Medium Compression Seamless Glove with Long Finger

**CcI III STRONG COMPRESSION GLOVES (34 -46 mmHg)**

CODE	PRODUCT
539D	Strong Compression Seamless Glove with Long Finger



**Seamless Glove**

## ■ SIZE TABLE


### Circumference of Hand(cm)

	SIZE						
Hand Circumference (cm)	1	2	3	4	5	6	7
cA	15-18	16,5-19	18-21	20-23	22-25	24-27	26-29
cB	15,5-19	17-20	18,5-22	20,5-24	22,5-26	24,5-28	26,5-30
cC	12-13	13-15	15-17	17-19	19-21	21-23	23-25

Length of Hand(cm)	
lAB	3-5
lAC	9-12

### Circumference of Fingers(cm)

		SIZE								
Finger Circumference (cm)		1	2	3	4	5	6	7	Length of Fingers	
Thumb	Top	5,5	6,0	6,5	7,0	7,5	8,0	9,0	Thumb	3,5-4,5
	Bottom	6,5	7,0	7,5	8,0	8,5	9,0	10,0		
Index Finger	Top	5,0	5,5	6,0	6,5	7,0	7,5	8,5	Index Finger	4,5-5,5
	Bottom	6,0	6,5	7,0	7,5	8,0	8,5	9,5		
Middle Finger	Top	5,0	5,5	6,0	6,5	7,0	7,5	8,5	Middle Finger	5,5-6,5
	Bottom	6,0	6,5	7,0	7,5	8,0	8,5	9,5		
Ring Finger	Top	4,5	5,0	5,5	6,0	6,5	7,0	8,0	Ring Finger	4,5-5,5
	Bottom	5,5	6,0	6,5	7,0	7,5	8,0	9,0		
Little Finger	Top	4,0	4,5	5,0	5,5	6,0	6,5	7,5	Little Finger	3,5-4,5
	Bottom	5,0	5,5	6,0	6,5	7,0	7,5	8,5		



**Choosing the size:** : To choose the exact size of Varimed® Lymphoedema Gloves, it is sufficient to determine the sizes of wrist (cC), hands(cA) and palms of hands elbow (cB). You can choose the suitable size from the measurements in the table. Do not use the glove in case of patient's whole arm and finger measurements do not correspond to sizes indicated in table.







ORTHOPAEDICS

# VARIMED ORTHOPAEDICS BACK SUPPORTS





## ■ VARIMED LUMBOSACRAL SUPPORT

### Indications:

- Lower backache (lumbago)
- Spondylolisthesis
- Spondylosis
- Degenerative spine diseases
- Spinal disc herniation
- Postoperative treatment after disc surgery
- Acute and chronic ache in the lumbosacral area

### Features:

- Four flexible steel splints at the back to provide extra support and stability to the lumbar area
- Anatomical design that fits both male and female
- Easy to put on and remove through the practical Velcro system
- Easy to put on by the hand straps on the belt
- Special flat knitted fabric for more elasticity, stability and comfortable to wear



**Lumbosacral Support**



**Lumbosacral Support With Belt**

## ■ CODES AND DESCRIPTIONS OF PRODUCTS

CODE	PRODUCT
<b>428</b>	Lumbosacral Support
<b>428B</b>	Lumbosacral Support with Belt

## ■ SIZE TABLE

	SIZE				
	S/1	M/2	L/3	XL/4	XXL/5
Waist Circumference(cm)	75-85	85-95	95-105	105-115	115-130

**Choosing the size:** To choose the exact size of Varimed® Elastic Corset, it is sufficient to determine the sizes of waist circumference. You can choose the suitable size from the measurements in the table.







## VARIMED ABDOMINAL CORSET



## ■ VARIMED Abdominal Basic 20 cm

### Indications:

- Abdominal stabilisation
- Support abdominal wall for healing after a surgery
- Abdominal muscular traumatism

### Features:

- Air permeable fabric
- Elastic fabric that avoids sweating
- Easy put and remove due to practical velcro system



**Abdominal Basic 20 cm**

## ■ VARIMED Abdominal Basic 26 cm

### Indications:

- Abdominal stabilisation
- Support abdominal wall for healing after a surgery
- Abdominal muscular traumatism

### Features:

- Air permeable fabric
- Elastic fabric that avoids sweating
- Easy put and remove due to practical velcro system



**Abdominal Basic 26 cm**

## ■ COLORS



**Beige**

## ■ SIZE TABLE

SIZE					
	S	M	L	XL	XXL
Waist Circumference (cm)	65-75	75-85	85-100	100-105	115-130

## ■ VARIMED Abdominal Plus 20 cm

### Indications:

- Abdominal stabilisation
- Support abdominal wall for healing after a surgery
- Abdominal muscular traumatism

### Features:

- Special flat knitted fabric for more elasticity, stability and comfortable to wear
- Easy put and remove due to practical velcro system

## ■ VARIMED Abdominal Plus 26 cm

### Indications:

- Abdominal stabilisation
- Support abdominal wall for healing after a surgery
- Abdominal muscular traumatism

### Features:

- Special flat knitted fabric for more elasticity, stability and comfortable to wear
- Easy put and remove due to practical velcro system

## ■ COLORS



Grey



## ■ SIZE TABLE

SIZE					
	S	M	L	XL	XXL
Waist Circumference (cm)	65-75	75-85	85-100	100-105	115-130



***Varimed Abdominal Basic***



***Varimed Abdominal Plus***



# VARIMED ORTHOPAEDICS KNEE SUPPORTS



## VARIMED LIGAMENT AND PATELLA SUPPORTED KNEE BRACE

### Indications:

- Osteoarthritis of knee joint
- Arthritis
- Patella instability
- Damage to meniscus
- Patella femoral pain syndrome
- Ligament instability

### Features:

- Silicone lined stays ensure nonslip grip for brace stability
- Spiral steel springs help support both side of knee
- Patella fixation support with silicone massaging pad for stabilisation, relief of the knee joint and rapid reduction of oedema
- Special flat knitted fabric for more elasticity, stability and comfortable to wear



## VARIMED LIGAMENT AND PATELLA SUPPORTED KNEE BRACE WITH OPEN PATELLA

### Indications:

- Patella chondropathy
- Meniscopathy
- Arthritis
- Meniscus Injuries
- Bursitis

### Features:

- Remove completely of the pressure on patella
- Silicone lined stays ensure nonslip grip for brace stability
- Spiral steel springs help support both side of knee
- Patella fixation support with silicone massaging pad for stabilization, relief of the knee joint and rapid reduction of oedema
- Special flat knitted fabric for more elasticity, stability and comfortable to wear



## ■ VARIMED PATELLA SUPPORTED KNEE BRACE

### Indications:

- Patella instability
- Damage to meniscus
- Patella femoral pain syndrome
- Patella chondropathy

### Features:

- Silicone lined stays ensure nonslip grip for brace stability
- Patella fixation support with silicone massaging pad for stabilisation, relief of the knee joint and rapid reduction of oedema
- Special flat knitted fabric for more elasticity, stability and comfortable to wear



## ■ VARIMED KNEE BRACE

### Indications:

- Provide resistance to a lateral blow when the knee is in full extension,
- Prevention of tendinous and articular accidents at work and practising sports,
- Rheumatic knee diseases,
- Bursitis
- Artritis

### Features:

- Special flat knitted fabric for more elasticity, stability and comfortable to wear
- Effective compression
- Anatomically formed



## ■ VARIMED PATELLA STRAP

### Indications:

- Patellar chondropathy
- Femoropatellar arthrosis
- Patellar Tendonitis

### Features:

- Soft bordering
- Pain relief
- No pressure in the hollow of the knee
- Soft patella pad supports patella tendon area



*Note: Standart Size*




## ■ CODES AND DESCRIPTIONS OF VARIMED KNEE SUPPORTS

CODE	PRODUCT
<b>423</b>	Knee Brace
<b>423D</b>	Patella Supported Knee Brace
<b>423A</b>	Ligament and Patella Supported Knee Brace with Open Patella
<b>423L</b>	Ligament and Patella Supported Knee Brace
<b>429P</b>	Patella Strap

## ■ SIZE TABLE

SIZE					
	XS	S	M	L	XL
Knee Circumference (cE) cm	30-34	34-38	38-42	42-47	47-53



**Choosing the Size:** To choose the exact size of Varimed® Knee Braces it is sufficient to determine the circumference of knee (cE). You can choose the suitable size from the measurements in the table.







# VARIMED ORTHOPAEDICS ANKLE SUPPORTS



## ■ VARIMED MALLEOLI SUPPORTED ANKLE BRACE

### Indications:

- Posttraumatic and postoperative irritations of ankle
- Ligament weakness
- Ligament Injuries
- Artroris
- Osteoarthritis
- Sprains
- Chronic instability

### Features:

- Two silicone pads that accelerate healing by massage effect
- Anatomically designed heel part supports foot comfort and ensures correct fit
- Special flat knitted fabric for more elasticity, stability and comfortable to wear



## ■ VARIMED LIGAMENT SUPPORTED ANKLE BRACE

### Indications:

- Protection of malleolus
- Ligament weakness
- Ligament injuries
- Chronic instability of ankle joint
- After sprains

### Features:

- Two plastic splints on the side of the ankle brace that guarding malleolus prevent sprains of introversion and extroversion of the ankle joint and minimizing its injury
- Special flat knitted fabric for more elasticity, stability and comfortable



## ■ VARIMED ACHILLES TENDON SUPPORTED ANKLE BRACE

### Indications:

- Acute and chronic inflammation of the Achilles tendon
- Achillodynia
- Achillobusitis
- Postoperative and post irritation

### Features:

- The compressive fabric and silicone pads that accelerate healing by massage effect
- Soft protective silicone pad for the Achilles tendon
- Special flat knitted fabric for more elasticity, stability and comfortable to wear



## ■ VARIMED ANKLE BRACE

### Indications:

- Distortions, bursitis, diarthrosis,
- Reduces the risk of incurring a future injuries,
- Prevents ankle sprains.

### Features:

- Special flat knitted fabric for more elasticity, stability and comfortable to wear
- Effective compression
- Anatomically formed



## ■ CODES AND DESCRIPTIONS OF VARIMED ANKLE SUPPORTS

CODE	PRODUCT
<b>421</b>	Ankle Brace
<b>421D</b>	Malleoli Supported Ankle Brace
<b>421L</b>	Ligament Supported Ankle Brace
<b>421T</b>	Achilles Tendon Supported Ankle Brace

SIZE					
	XS	S	M	L	XL
Ankle Circumference (cB) cm	18-20	20-23	23-26	26-29	29-32



**Choosing the Size:** To choose the exact size of Varimed® Elastic Ankle Braces it is sufficient to determine the circumference of ankle (cB). You can choose the suitable size from the measurements in the table.





## VARIMED ORTHOPAEDICS ARM SUPPORT



## ■ VARIMED SHOULDER BRACE

### Indications:

- Shoulder Stabilisation
- Osteoarthritis
- Chronic, postoperative and posttraumatic irritations of the soft tissues in the shoulder area.
- Humerus fractures

### Features:


- Active support for treatment of the shoulder
- Anatomically design covers the upper arm and shoulder
- Special flat knitted fabric for more elasticity, stability and comfortable to wear
- The simple strap system makes the shoulder support quick and easy to put on



## ■ CODES AND DESCRIPTIONS OF VARIMED SHOULDER BRACE

CODE	PRODUCT
<b>4240</b>	Shoulder Brace

SIZE					
	S	M	L	XL	XXL
Upper part of arm Circumference (cm)	23-25	26-28	29-31	32-34	35-38



# VARIMED ORTHOPAEDICS ELBOW SUPPORTS



## ■ VARIMED ELBOW BRACE WITH SILICONE PADS

### Indications:

- Chronic, postoperative and posttraumatic irritations of the elbow
- Epycondylitis
- Dislocation of the elbow
- Tendomyopathy

### Features:

- Anatomic design to fit the elbow area
- Two silicone pad that accelerate healing by massage effect
- Special flat knitted fabric for more elasticity, stability and comfortable to wear



## ■ VARIMED ELBOW BRACE

### Indications:

- Chronic, postoperative and posttraumatic irritations of the elbow
- Articular effusion
- Artrosis
- Sprain
- Contusion

### Features:

- Special flat knitted fabric for more elasticity, stability and comfortable to wear
- Effective compression
- Anatomically formed



## ■ VARIMED EPICONDYLITIS BANDAGE

### Indications:

- Epicondylitis
- Chronic or relapsing epicondylopathy

### Features:

- Elastic band that allows pressure adjustment
- One universal size
- The same design for right and left elbow
- Padded inside




Note: Standard Size

## ■ CODES AND DESCRIPTIONS OF VARIMED ELBOW SUPPORTS

CODE	PRODUCT
424	Elbow Brace
424D	Elbow Brace with Silicone Pads
427	Epicondylitis Bandage

## ■ SIZE TABLE

SIZE					
	XS	S	M	L	XL
5 cm. Lower part of the elbow (cD) cm	20-23	23-26	26-29	29-32	32-35



**Choosing the Size:** To choose the exact size of Varimed® Elbow Braces it is sufficient to determine the circumference 5 cm lower part on elbow (cD). You can choose the suitable size from the measurements in the table.





# VARIMED ORTHOPAEDICS HAND SUPPORTS



## ■ VARIMED THUMB SPLINT

### Indications:

- Chronic, postoperative and post traumatic irritations
- Degenerative conditions of the thumb joint
- Arthrosis of the carpometacarpal joint of the thumb
- Capsule and ligament lesions of the carpometacarpal joint

### Features:

- Immobilisation and stabilisation of the thumb due to formable aluminium splint
- Produced mainly by skin friendly and breathable fabric
- Standart size



## ■ VARIMED WRIST BRACE WITH SPLINT

### Indications:

- Chronic, postoperative and post traumatic irritations
- Carpal arthroris
- Carpal tunel syndrome
- Tenosynovitis
- Wrist instability

### Features:

- Immobilisation and stabilisation of the wrist due to formable aluminium splint
- Anatomic design to fit the wrist
- Produced mainly by skin friendly and breathable fabric
- Fit for both left and right hand



## ■ VARIMED WRIST SPLINT WITH THUMB GRIP

### Indications:

- Chronic, postoperative and posttraumatic irritations
- Severe distortion of the wrist and thumb
- Postoperative fixation of the wrist and thumb
- Degeneration in the wrist and thumb area

### Features:

- Immobilisation and stabilisation of the wrist and thumb due to formable aluminium splint
- Anatomic design to fit the wrist and thumb
- Produced mainly by skin friendly and breathable fabric




## ■ CODES AND DESCRIPTIONS OF VARIMED HAND SUPPORTS

CODE	PRODUCT
<b>426</b>	Wrist Brace with Splint
<b>426B</b>	Thumb Splint (Standart Size)
<b>426BR</b>	Wrist Splint with Thumb Grip (Right Hand)
<b>426BL</b>	Wrist Splint with Thumb Grip (Left Hand)

## ■ SIZE TABLE

SIZE					
	XS	S	M	L	XL
Wirst Circumference (cm)	13-15	15-17	17-19	19-21	21-23



**Choosing the Size:** To choose the exact size of Varimed® Elastic Elbow Braces it is sufficient to determine the wrist circumference (cC). You can choose the suitable size from the measurements in the table.





## VARIMED ELASTIC BRACES



## ■ SPORTS INJURIES AND ELASTIC BRACES

The principles of sport medicine can be applied to the treatment of many musculoskeletal injuries which resemble sports injuries but have different causes. For example tennis elbow can be caused by carrying a suitcase, turning a screw or opening a stuck door and runner's knee can be caused by excessive inward rolling of the foot (pronation) during walking.

Sports injuries are caused by incorrect training methods, structural abnormalities that stress certain parts of the body more than others and weakness of muscles, tendons and ligaments. Many of these injuries are caused by chronic wear and tear which results from repetitive motion stressing susceptible tissue.

The most common cause of muscle and joint injuries are incorrect training methods: The exerciser doesn't allow for adequate recovery after a workout or doesn't stop exercising when pain develops.

Structural abnormalities can cause a person to be susceptible to sports injuries by stressing parts of the body unevenly. For example when the legs are unequal in length, greater force is exerted on the hip and knee of the longer leg. In general, the pain disappears when the activity is stopped but reoccurs each time the same exercise intensity is reached.

Muscles, tendons and ligaments tear when exposed to forces greater than their inherent strength. They may be injured if they are too weak or too tight for the exercise being attempted.

Some authors have determined favorable effects of compression stockings, tights or sleeves on venous hemodynamics, arterial perfusion, deeper tissue oxygenation and muscle oscillation with corresponding effects on lactate kinetics and muscular soreness.

Prophylactic ankle braces generally reduce the risk of injury during sports and are used in treatment of acute ankle sprain.

Many biomechanical studies have shown that prophylactic knee braces can provide 20% to 30% greater resistance to a lateral blow with the possibility that the anterior cruciate ligament is given even greater protection than the medial collateral ligament, additionally the elastic support helps to limit swelling.

Other scientific studies show that elastic calf braces may provide positive changes in performance of repeated sprints and maximal isokinetic exercise. Also they may reduce pain perception. Elastic braces with dynamic gradient compression accelerate the flow of oxygen-rich blood to the muscles to increase power and endurance and the removal of lactic acid from the muscles. This leads to an increase of stability and improves the management of body temperature and moisture.

## ■ VARIMED ELASTIC ANKLE BRACE

### Indications:

- Distortions, bursitis, diarthrosis,
- Reduces the risk of incurring a future injuries,
- Helps to prevent ankle sprains

### Features:

- Effective compression
- Anatomically formed
- Highly elastic



## ■ VARIMED ELASTIC CALF BRACE

### Indications:

- Tendon and muscle tensions,
- For torn muscles, post-traumatic rehabilitation,
- Prevent DVT prophylaxis at the calf,

### Features:

- Effective compression
- Anatomically formed
- Highly elastic



## ■ VARIMED ELASTIC KNEE BRACE

### Indications:

- Provide resistance to a lateral blow when the knee is in full extension,
- Prevention of tendinous and articular accidents at work and practising sports,
- Rheumatic knee diseases,

### Features:

- Effective compression
- Anatomically formed
- Highly elastic



## ■ VARIMED ELASTIC ELBOW BRACE

### Indications:

- Chronic, postoperative and posttraumatic irritations of the elbow
- Articular effusion
- Arthrosis
- Sprain
- Contusion

### Features:

- Effective compression
- Anatomically formed
- Highly elastic



## ■ CODES AND DESCRIPTIONS FOR VARIMED ELASTIC BRACES

CODE	PRODUCT
<b>411</b>	Elastic Ankle Brace
<b>412</b>	Elastic Calf Brace
<b>413</b>	Elastic Knee Brace
<b>414</b>	Elastic Elbow Brace

## ■ SIZE TABLE

### Elastic Ankle Brace

SIZE					
	XS	S	M	L	XL
Ankle Circumference (cB) cm	18-20	20-23	23-26	26-29	29-32

### Elastic Calf Brace

SIZE					
	XS	S	M	L	XL
Calf Circumference (cC) cm	28-32	32-35	35-38	38-41	41-45

### Elastic Knee Brace


SIZE					
	XS	S	M	L	XL
Knee Circumference (cE) cm	30-34	34-38	38-42	42-47	47-53



## Elastic Elbow Brace

SIZE					
	XS	S	M	L	XL
5 cm.Lower part of the elbow (cD) cm	20-23	23-26	26-29	29-32	32-35

**Choosing the Size:** To choose the exact size of Varimed® Elastic Braces it is sufficient to determine the circumference of ankle (cB) for ankle braces, circumference of calf (cC) for calf braces, circumference of knee (cE) for knee braces and circumference 5 cm lower part on elbow (cD) for elbow braces. You can choose the suitable size from the measurements in the table.



## VARIMED Dynamic SUPPORTS



## ■ VARIMED Dynamic SPORT STOCKING

Varimed Dynamic Sports Stockings have a graduated profile with maximum compression at ankle, decreasing towards the thigh, thus improve the oxygen supply of muscles and ensures more energy. The inner surface of the sock that touches the skin is 100% cotton.

### Features:

- Sporty look
- Anatomic fit
- Microfilament yarns give the stocking excellent breathability
- High moisture transportation for pleasant wearing comfort
- Perfect fit because of soft stretch material,
- Ideal for all running sports, ball games and athletics,
- For both men and women

### Benefits:

- Anatomical support for higher performance
- Reduced feeling of exhaustion,
- Compression pressure increases blood circulation and venous performance in motion.



## ■ VARIMED Dynamic BACK SUPPORT

### Indications:

- Light Support for muscular stabilisation of the lumbar
- Lower backache (lumbago)

### Features:

- Flat knit
- Anatomical design that fits both male and female
- Easy to put on and remove through the practical Velcro system
- Special flat knitted fabric for more elasticity, stability and comfortable to wear
- At work, at home, in sports



## ■ VARIMED Dynamic KNEE SPORT

### Indications:

- Provide resistance to a lateral blow when the knee is in full extension,
- Prevention of tendinous and articular accidents at work or practising sports,
- Rheumatic knee diseases,

### Features:

- Special flat knitted fabric for more elasticity, stability and comfortable to wear
- Effective compression
- Anatomically formed



## ■ VARIMED Dynamic ANKLE SUPPORT

### Indications:

- Distortions, bursitis, diarthrosis,
- Reduces the risk of incurring a future injuries,
- Prevents ankle sprains

### Features:

- Special flat knitted fabric for more elasticity, stability and comfortable to wear
- Effective compression
- Anatomically formed



## ■ VARIMED Dynamic ELBOW SUPPORT

### Indications:

- Chronic, postoperative and posttraumatic irritations of the elbow
- Articular effusion
- Artrosis
- Sprain
- Contusion

### Features:

- Special flat knitted fabric for more elasticity, stability and comfortable to wear
- Effective compression
- Anatomically formed





## ■ VARIMED Dynamic HAND/WRIST SUPPORT

### Indications:

- Prevention of overstraining and sprains
- Wrist instability

### Features:

- Anatomic design to fit the wrist
- Special flat knitted fabric for more elasticity, stability and comfortable to wear
- Fit for both left and right hand



*Hand/Wrist Support*

## ■ VARIMED Dynamic WRIST SUPPORT

### Indications:

- Excellent support effect
- Helps intensify the sense of joint stability

### Features:

- Anatomic design to fit the wrist
- Special flat knitted fabric for more elasticity, stability and comfortable to wear
- Fit for both left and right hand



*Wrist Support*

## ■ CODES AND DESCRIPTIONS FOR VARIMED Dynamic SUPPORTS

CODE	PRODUCT
<b>428C</b>	Varimed Dynamic Back Support
<b>423C</b>	Varimed Dynamic Knee Support
<b>421C</b>	Varimed Dynamic Ankle Support
<b>424C</b>	Varimed Dynamic Elbow Support
<b>426C</b>	Varimed Dynamic Hand/Wrist Support
<b>429C</b>	Varimed Dynamic Wrist Support
<b>91P1</b>	Varimed Dynamic Sport Stocking

## ■ SIZE TABLE

### Varimed Dynamic Ankle Brace Support

SIZE				
	S	M	L	XL
Ankle Circumference (cB) cm	19-22	22-25	25-28	28-32

### Varimed Dynamic Knee Brace Support

SIZE				
	S	M	L	XL
Knee Circumference (cE) cm	32-37	37-42	42-47	47-52

### Varimed Dynamic Elbow Support

SIZE				
	S	M	L	XL
5 cm.Lower part of the elbow (cD) cm	20-23	23-27	27-31	31-35

### Varimed Dynamic Hand / Wrist Support

SIZE				
	S	M	L	XL
Wrist Circumference (cC) cm	13-15	15-18	18-21	21-24

### Varimed Dynamic Wrist Splint Support


SIZE				
	S	M	L	XL
Wrist Circumference (cC) cm	13-15	15-18	18-21	21-24

### Varimed Dynamic Back Support

SIZE				
	S	M	L	XL
Waist Circumference (cm)	65-80	80-95	95-110	110-125

### Varimed Dynamic Sport Stocking

	SIZE				
Leg Circumference (cm)	1	2	3	4	5
cB	18-20	20-23	23-26	26-29	29-32
cC	28-34	31-37	34-41	37-44	41-48



**Choosing the Size:** To choose the exact size of Varimed Dynamic Supports it is sufficient to determine the circumference of ankle (cB) for ankle braces, circumference of knee (cE) for knee braces, circumference 5 cm lower part on elbow (cD), circumference of ankle (cB) and calf (cC) for Sport Stockings, circumference of wrist (cC) for Hand/Wrist Supports and circumference of waist for Back Supports. You can choose the suitable size from the measurements in the table.



## VARIMED Thermo SUPPORTS





## ■ VARİMED Thermo ELASTIC WAIST CORSET

### Indications:

- To keep the waist, abdomen and kidneys warm
- Spasm of muscles
- Pains due to rheumatism and sciatica

### Features:

- Excellence protection from cold in the lumbar spinal area



## ■ VARİMED Thermo ELASTIC KNEE BRACE

### Indications:

- Rheumatism
- Sciatica
- Any kinds of knee pains

### Features:

- Excellence protection from cold in the knee area



## ■ CODES AND DESCRIPTIONS FOR VARİMED Thermo SUPPORTS

KOD	ÜRÜN
428Y	Thermo Elastic Waist Corset
428YC	Thermo Elastic Adjustable Waist Corset
423Y	Thermo Elastic Knee Brace

## ■ SIZE TABLE

### Thermo Elastic Knee Brase

SIZE					
	S	M	L	XL	XXL
Knee Circumference(cE) cm	30-34	34-38	38-42	42-47	47-53

### Thermo Elastic Waist Corset

SIZE					
	S	M	L	XL	XXL
Waist Circumference (cm)	75-85	85-95	95-105	105-115	115-130





## ACCESSORIES



## **VARIMED DONNING AID FOR MEDICAL COMPRESSION STOCKINGS**

Quick and simple donning of all type medical compression Stockings



## **VARIMED DRESSING AID FOR MEDICAL COMPRESSION STOCKINGS - FABRIC**

Ideal donning aid for open toe medical compression stockings



## **VARIMED DRESSING AID FOR MEDICAL COMPRESSION ARMSLEEVE FABRIC**

Ideal donning aid for medical compression armsleeves





## ■ SCIENTIFIC EVIDENCE OF COMPRESSION TREATMENT

In the following the results of two scientific studies that show the effects of compression treatment on chronic vascular diseases:

\*This consensus study as per editor Hugo Partschi is the result of an international expert meeting organized by the International Union of Phlebology (IUP) in August 2003 (Table 1). Recommended levels of compression in treatment of A, B and C are rated at three levels. The higher statistic significant level, large randomize control trials and metaanalyses with homogenous results were assigned to a recommendation level A, the higher statistic significant level, smaller or the higher statistic significant level, single randomize control trials to a level B and observational studies or the opinion of the experts were taken as level C. The following table summarizes the clinical indications for compression class (A, I, II, III) according to the European CEN regulation of recommendation for stockings. Chronic venous disorders are differentiated following the CEAP classification: C0 no visible varices, C1 teleangectasias, C2 large varicose veins, C3 oedema, C4 skin changes, C5 ulcer scar, C6 venous ulcer.

Table 1. Scientific Evidence Of Compression Treatment, Hugo Partschi (2003)

INDICATIONS	Ccl A	Ccl A	Ccl II	Ccl III
C0, C1 symptomatic	B	B		
C1			B	B
C2 asymptomatic			C	
C2 symptomatic				C
C2 pregnancy			B	
C2 After surgery			C	C
C2 After sclerotherapy				B-C
C3			B	
C4b (LDS)			B	
C5			B	B
C6			B	
DVT prevention	A-B	A-B		
Flights	B	B		
DVT therapy			B	B
PTS prevention			B	A
Lymphoedema			C	C

DVT= Deep Vein Thrombosis, PTS= Postthrombotic Syndrome

\*\* Marcondes Figuieredo<sup>2</sup>, summarized the scientific evidence in the meeting of International Compression Club in Paris in November 2007 (Table 2).

Table 2- Scientific Evidence Of Compression Treatment, Marcondes Figueiredo (2009)

INDICATIONS	Compression (mm Hg)		
	10-20	20-30	30-40
C0 s	1B		
C1 s	1B		
C1 s post-sclerotherapy		1B	
C2 a, s	1B	2B	
C2 s pregnancy	1B	1B	
C3 prevention	1B		
C3 therapy		2B	
C4 B		1B	1B
C5		2B	1A
C6			1B
VTE Prevention	1A		
VTE Therapy			1B
PTS Prevention			1A
PTS Therapy			2B
Lymphoedema Therapy			1B

a=asymptomatic; s=symptomatic; PTS=post-thrombotic syndrome; VTE=venous thromboembolism

# CERTIFICATE

## The company

**Varimed Medikal San. Ve Tic. Ltd. Sti.**  
Terazidere Mah. Davutpasa Ciftehavuzlar Cad. No. 28-30  
34035 Bayrampasa - Istanbul, TURKEY

is granted authorisation according to STANDARD 100 by OEKO-TEX® to use  
the STANDARD 100 by OEKO-TEX® mark, based on our test report  
**18.0.45019**



## for the following articles:

Knee length and thigh length anti-embolism stockings, leg ulcer kits, support stockings, medical compression stockings, lymphoedema gloves, lymphoedema armsleeves and lymphoedema stockings, diabetic socks, elastic braces, sport stockings, made of yarn dyed knitted fabrics (partly with flat band) produced from polyamide, polyamide/elastane, cotton/elastane, polyester, polyester/elastane and metalized yarn and mixtures among each other as well as white and dyed woven fabric produced from polyester, partly finished with biological active products accepted by OEKO-TEX®, incl. accessories (hook and loop fasteners, woven ribbons with silicone coating, sewing threads, printed labels [black on white], elastic yarns, zippers); produced by using materials certified according to STANDARD 100 by OEKO-TEX®.

The results of the inspection made according to STANDARD 100 by OEKO-TEX®, Appendix 4, **product class II** have shown that the above mentioned goods meet the human-ecological requirements of the STANDARD 100 by OEKO-TEX® presently established in Appendix 4 for products with direct contact to skin.

The certified articles fulfil requirements of Annex XVII of REACH (incl. the use of azo colourants, nickel release, etc.), the American requirement regarding total content of lead in children's articles (CPSIA; with the exception of accessories made from glass) and of the Chinese standard GB 18401:2010 (labelling requirements were not verified).

The holder of the certificate, who has issued a conformity declaration according to ISO 17050-1, is under an obligation to use the STANDARD 100 by OEKO-TEX® mark only in conjunction with products that conform with the sample initially tested. The conformity is verified by audits.

**The certificate 14.HTR.36217 is valid until 31.05.2019**

Boennigheim, 06.08.2018

  
Dipl.-Ing. (FH) Elisabeth Panian  
Head of Certification Body OEKO-TEX®







2018

## CERTIFICATE OF REGISTRATION

*This certifies that:*

**VARIMED MEDIKAL SAN. VE TIC. LTD. STI**  
**Terazidere Mah.**  
**Davutpasa Ciftehavuzlar Cad. No. 28 - 30**  
**Bayrampasa Istanbul, TURKEY TR-34035**

is registered with the U.S. Food and Drug Administration for FY 2018 pursuant to Title 21, 807 et seq. of the United States Code of Federal Regulations:

Establishment Registration:	<b>3010375113</b>
DUNS No.:	<b>36-564-7414</b>
Device Classification Name:	<b>STOCKING, MEDICAL SUPPORT (TO PREVENT POOLING OF BLOOD IN LEGS)</b>
Product Code:	<b>DWL</b>
Regulation Number:	<b>880.5780</b>
Official Correspondent and U.S. Agent:	<b>Registrar Corp</b> 144 Research Drive, Hampton, Virginia, 23666, USA Telephone: +1-757-224-0177 • Fax: +1-757-224-0179


*Registrar Corp will confirm that such registration remains effective upon request and presentation of this certificate until the end of the year stated above, unless said registration is terminated after issuance of this certificate. Registrar Corp makes no other representations or warranties, nor does this certificate make any representations or warranties to any person or entity other than the named certificate holder, for whose sole benefit it is issued. This certificate does not denote endorsement or approval of the certificate-holder's device or establishment by the U.S. Food and Drug Administration. Registrar Corp assumes no liability to any person or entity in connection with the foregoing.*

*Pursuant to 21 CFR 807.39, "Registration of a device establishment or assignment of a registration number does not in any way denote approval of the establishment or its products. Any representation that creates an impression of official approval because of registration or possession of a registration number is misleading and constitutes misbranding."*

*The U.S. Food and Drug Administration does not issue a certificate of registration, nor does the U.S. Food and Drug Administration recognize a certificate of registration. Registrar Corp is not affiliated with the U.S. Food and Drug Administration.*

**Registrar Corp** ★

144 Research Drive, Hampton, Virginia, 23666, USA  
Telephone: +1-757-224-0177 • Fax: +1-757-224-0179  
info@registrarcorp.com • www.registrarcorp.com

  
Russell K. Statman  
Executive Director  
Registrar Corp

Dated: December 20, 2017

# Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016 & EN ISO 13485:2016

This is to certify that:

Varimed Medikal Sanayi  
Ve Ticaret Ltd. Sti.  
Terazidere Mah.  
Çiftelavuzlar Cad. No:28  
Bayrampasa  
İstanbul  
34035  
Turkey

Holds Certificate Number:

**MD 563572**

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 & EN ISO 13485:2016 for the following scope:

The design, manufacture and sales of non-sterile medical compression stockings, anti-embolism stockings, diabetic socks, lymphoedema armsleeves, elastic ankle braces, elastic calf braces, elastic knee braces, leg ulcer kits, travel stockings and knitting orthoses and the import and selling of medical gloves.



For and on behalf of BSI:

Stewart Brain, Head of Compliance & Risk - Medical Devices

Original Registration Date: 2010-08-20

Latest Revision Date: 2019-02-08

Effective Date: 2016-08-20

Expiry Date: 2019-08-19

Page: 1 of 1



...making excellence a habit™

This certificate was issued electronically and remains the property of BSI and is bound by the conditions of contract.

An electronic certificate can be authenticated [online](https://www.bsigroup.com/ClientDirectory).

Printed copies can be validated at [www.bsigroup.com/ClientDirectory](https://www.bsigroup.com/ClientDirectory)

Information and Contact: BSI, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes MK5 8PR Tel: +44 345 080 9000  
BSI Assurance UK Limited, registered in England under number 7805321 at 389 Chiswick High Road, London W4 4AL, UK.  
A Member of the BSI Group of Companies.





## DECLARATION OF CONFORMITY

**MANUFACTURER:**

Name: Varimed Medikal Sanayi ve Ticaret Limited Şirketi

Address: Terazidere Mah. Davutpaşa Çiftelavuzlar Cad. No:28-30 Bayrampaşa 34035 İSTANBUL/TÜRKİYE

Phone:+90 212 493 02 35 (pbx) • Fax: +90 212 544 12 14

Web: www.varimed.com.tr • mail: info@varimed.com.tr

We, the undersigned, hereby declare that the described products meet the essential requirements listed in Medical Device Directive 93/42/EEC as ammended by 2007/47/EC.

This declaration is based on:

Essential Requirements Listed in Directive 93/42/EEC

Approval of a QMS System:

ISO 13485:2016

This declaration issupported by:

Varimed Medikal San. ve Tic. Ltd. Şti.-Technical File

Products:

Medical Compression Stockings,

Anti-EmbolismStockings,

Lymphoedema Products,

Elastic Support Products,

Leg Ulcer Kit,

Support Stockings

Travel Stockings,

Diabetic Socks,

Classification: Class 1, Rule 1

**SIGNED ON BEHALF OF THE MANUFACTURER**

Name : Dr. M. VOLKAN ÖZTÜRKATALAY

Position : GENERAL MANAGER

Place : İSTANBUL-TÜRKİYE

Date : 12.02.2019

Signature :

## VARIMED Production Line



## ■ VARIMED PRODUCTION LINE



### **Raw Material**

The best quality raw materials are supplied for each product group and their quality control is carried out in detail.



### **Knitting Process**

Products are knitted in a circular knitting or flat knitting machine according to the property.







### **Sewing Process**

All sewing operations such as topband stitching, net stitching, edge stitching and label stitching are performed.



### **Quality Control Process**

The products are subjected to pressure and elasticity tests according to their properties and are subjected to a detailed control against possible production defects.



### **Packing & Shipping**

The products are placed in high quality and Varimed designed boxes with detailed user guides. All the descriptive information of the products such as size, compression and color are attached to the box by a label and delivered to our customers.







**VARİMED MEDİKAL SAN. VE TİC. LTD. ŞTİ.**

Terazidere Mah. Davutpaşa Çiftehavuzlar Cad. No:28-30 Bayrampaşa 34035 İstanbul - Türkiye

Phone: +90 212 493 02 35 • Fax: +90 212 544 12 14

**info@varimed.com.tr • www.varimed.com.tr**



varimedmedikal



varimedmedikal



varimedmedikal